

# Timers

## digital timers

### Timer, number keys, Fisherbrand®, Traceable®

Fisherbrand

- Individually numbered keys make setting time faster than units with up/down buttons
- Counts down with alarm and counts up with time-out feature
- 4 digit LCD display shows the time remaining in minutes and seconds
- Stopwatch mode: counts up to 59min, 59s and then repeats
- Countdown mode: 1s increments from 99min, 59s
- Alarm: sounds at zero for 60s
- Big digit, 9.5mm high display
- Magnetic back, flip open easel for standing, a spring fastener for clipping to a lab coat, and a hole for a lanyard
- Accuracy: 0.01%
- Size: 63mm x 57mm x 13mm
- Mass: 42g
- Supplied with a 2 year, continuous use, silver oxide battery



Catalogue No	Legacy No	Description
11785863	-	Alarm/stopwatch timer with number keys

### Timer, LCD, nano, count down/up, Fisherbrand®, Traceable®

Fisherbrand

- Counts down/up, to/from 24hr with alarm
- Two timing channels
- 10mm and 5mm digits
- Memory
- Clock
- Back clip with magnet
- Dimensions: 60mm x 50mm
- Battery included



Catalogue No	Legacy No	Description
12695296	-	Timer, counts up/down, with alarm

#### Spare part

Catalogue No	Legacy No	Description
12727565	-	Button battery, 1.5V, 100mAh (LR44)

### Timer, flashing LED alert, big digit, Fisherbrand®, Traceable®

N

Fisherbrand

- Massive 5cm high display
- 3 countdown times can be stored to memory
- Counts down from 99min 59s to 1s, green light flashes 10s from zero
- When zero is reached, red light flashes, alarm sounds and timer starts counting up
- Alarm sets at loud, low or off (light flashing only)
- Magnetic back, flip open easel, spring fastener
- Accuracy: 0.01%
- Size: 88.9mm x 57.2mm x 15.8mm
- Mass: 71g
- Supplied with AAA alkaline batteries



Catalogue No	Legacy No	Description
12922830	-	Flashing LED alert big digit alarm timer

### Timer, 100 hour, Fisherbrand®, Traceable®

Fisherbrand

- Times up/down to/from 99hr, 59min, 59s
- Channels are four countdown, one clock and one stopwatch
- Time-out mode
- Alarm sounds for 1 min, display flashes 'TIME IS UP'
- Accuracy: 0.01%
- Mini-size: 70mm x 60mm x 16mm
- Supplied with batteries (expected battery life 2yr)



Catalogue No	Legacy No	Description
11512793	-	100 hour timer

#### Spare part

Catalogue No	Legacy No	Description	Pack qty
12326968	-	Batteries, alkaline, 130mAh, 1.5V (LR44)	4

For further Traceable® timers please see the Fisher Scientific Traceable® catalogue or visit our website