



Tripods, triangular and circular top



Zinc plated steel with 210mm height, 6mm diameter legs.

Catalogue No	Legacy No	Alt. No	Description
11779806	-	3001	Triangular, 125mm side length
12966916	-	3002	Triangular, 150mm side length
12956916	-	3004	Circular, 150mm O.D.



Trolleys, three tier, polypropylene

Manufactured from polypropylene and available in three different configurations, this range of trolleys is ideal for use in the laboratory.

Catalogue No	Legacy No	Alt. No	Dimensions [w x d x h], mm	Description	Colour
13110183	-	6420/M/BI	460 x 320 x 830	3 shelves and 1 tray	White
13120183	-	6420/M/GR	460 x 320 x 830	3 shelves and 1 tray	Grey
13130183	-	6425/M/BI	460 x 320 x 810	3 shelves	White
13140183	-	6425/M/GR	460 x 320 x 810	3 shelves	Grey
13150183	-	6410/M/BI	460 x 320 x 830	3 shelves, doors and 1 tray	White
13160183	-	6410/M/GR	460 x 320 x 830	3 shelves, doors and 1 tray	Grey

Spare parts

Catalogue No	Legacy No	Alt. No	Dimensions [w x d x h], mm	Description	Colour
13170183	-	6430/M/BI	460 x 320 x 270	Spare shelf	White
13180183	-	6430/M/GR	460 x 320 x 270	Spare shelf	Grey



Trolleys, mobile laboratory, Bobby®

Multi-purpose, mobile storage unit. This versatile unit has a unique combination of pivoting trays, shelves and handy recesses. Manufactured from white ABS plastic (5mm thick). The 2/2 models are low level, 2 tier units with 2 pivoting trays, the 4/4 and 4/6 models are high level 4 tier units, with either 4 or 6 pivoting trays. The top work surface of all models measures 410mm x 430mm. Mounted on heavy duty rubber castors giving total stability.

Catalogue No	Legacy No	Description	Height, mm
12375699	-	Trolley 2/2	530
13120193	-	Trolley 4/4	930
13130193	-	Trolley 4/6	930



Trolley

With epoxy powder coated metal frame and 76mm castors capable of transporting loads up to 100kg (evenly distributed). The 80mm deep polyethylene trays are removable. The trolley can be dismantled for transport or storage.

Catalogue No	Legacy No	Dimensions [w x d x h], mm
11346914	-	810 x 470 x 1,050 (including handles)

Spare part

Catalogue No	Legacy No	Description
11350685	-	Tray

For mini-dolly folding trolleys please see the Traceable® Products Catalogue.