

# Tubes

## test tubes



### Caps, phenolic, rubber liner, for screw thread test tubes, Fisherbrand®



**Fisherbrand®**

Made of a special phenolic material (wood filled) and very resistant to the effects of temperature and steam encountered in autoclaving. A cemented-in rubber liner is provided, and the cement used has been selected to retain its adhesive properties during autoclaving.

Catalogue No	Legacy No	GPI thread finish	Pack qty
11547133	-	13-415	1,000
12337279	-	15-415	1,000
11567133	-	18-415	1,000



### Caps, polypropylene, linerless, for screw thread test tubes, Fisherbrand®



**Fisherbrand®**

One-piece construction with a unique inner sealing ring.

Catalogue No	Legacy No	GPI thread finish	Pack qty
11577133	-	13-415	1,000
11587133	-	15-415	1,000
11597133	-	18-415	500



### Caps, polypropylene, welded PTFE/silicone liner, for screw thread test tubes, Fisherbrand®



**Fisherbrand®**

Ideal for repeated autoclaving.

The welded liner eliminates glue contamination.

Catalogue No	Legacy No	GPI thread finish	Pack qty
11522672	-	13-415	288
11532672	-	15-415	288
11542672	-	18-415	288



### Test tubes, polystyrene, round base, Thermo Scientific

**Thermo**  
SCIENTIFIC

Transparent polystyrene, round base, suitable for samples from 0.4mL to 4.9mL.

Catalogue No	Legacy No	Alt. No	Dimensions, mm	Capacity, mL	Pack qty
11748868	-	RT15	40mm x 6mm	0.4	10,000
11738858	-	112	40mm x 11mm	2.3	1,000
13015723	-	RT20	50mm x 6mm	1.0	10,000
11708988	-	RT25	65mm x 10mm	2.7	6,000
11778908	-	RT30	75mm x 12mm	4.9	3,500

#### Caps

Catalogue No	Legacy No	Alt. No	Description	Pack qty
11768908	-	C2A	For 11778908, 4.9mL capacity tube	3,500