

## Vials, storage, with PTFE lined caps, Vial File™



Open top caps are lined with PTFE faced silicone rubber. Stored in Vial File™.



Catalogue No	Legacy No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
11571205	-	224950	1.5	12	35	8	425	60
10206112	-	224952	4	15	48	13	425	40

## Vials, V, screw top with PTFE lined solid top caps



Solid top black phenolic screw cap with PTFE faced 14B white rubber liner.



Catalogue No	Legacy No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
11506305	-	W986254NG	1.0	13	44	13	425	12
11596295	-	W986253NG	0.3	13	35	13	425	12
11516305	-	W986257NG	3.0	20	50	20	400	12

## Vials, V, screw top, graduated, with PTFE lined open top caps



Graduated vial with open top black phenolic screw cap and PTFE faced silicone liner.



Catalogue No	Legacy No	Alt. No	Capacity, mL	Diameter, mm	Height, mm	Closure size, mm	GPI thread finish	Pack qty
11556305	-	986293	0.3	13	32	13	425	12
11566305	-	986294	1.0	13	41	13	425	12
13111813	-	986298	2.0	20	41	20	400	12
11576305	-	986297	3.0	20	47	20	400	12
10299974	-	986299	5.0	20	62	20	400	12

For chromatography vials and closures, please see pages 318 to 342

