

Vials, HDPE, liquid scintillation, with caps, 7mL, Fisherbrand®



Fisherbrand®

7mL smaller capacity vials require less solvent. HDPE provides low background counts and excellent light transmission, resists swelling caused by solvent degradation. White caps fit GPI 15-425 thread finish and are suitable for marking.

Packs of 1,000 vials are shrinkwrapped in four trays of 250, with caps separately packed in four polyethylene bags of 250.

Packs of 2,000 vials are bulk packed in one polyethylene bag, with caps packed in two polyethylene bags of 1,000.

Catalogue No	Legacy No	Cap material	Cap style	Packaging type	Pack qty
12371599	-	Polyethylene	No liner	Caps separate	1,000
12644347	-	Polyethylene	No liner	Caps separate	2,000



Vials, HDPE, liquid scintillation, with caps, 20mL, Fisherbrand®



Fisherbrand®

20mL vials fit all common counters.

High density polyethylene (HDPE) provides low background counts, excellent light transmission and resists leaking and swelling caused by solvent degradation. White caps fit GPI 22-400 thread, and are suitable for marking. Vials with caps attached are shrinkwrapped in five trays of 100 per case. Where caps are separate, vials are shrinkwrapped in five trays of 100 per case, and caps are packed in five polybags of 100 each (except 12361599, which is bulk packed with 1,000 in one polybag, with caps in a separate polybag of 1,000)

Vial dimensions, with cap, (O.D. x h), mm: 28 x 61.

Catalogue No	Legacy No	Cap material	Cap style	Packaging type	Pack qty
12321599	-	Urea	Cork backed metal foil liner	Caps separate	500
12365573	-	Urea	Cone shaped plastic liner	Caps separate	500
12341599	-	Polypropylene	Pulp backed metal foil liner	Caps separate	500
12351599	-	Polypropylene	Pulp backed metal foil liner	Caps attached	500
12361599	-	Polypropylene	Pulp backed metal foil liner	Bulk pack	1,000

Vials, HDPE, liquid scintillation


WHEATON

Wheaton polyethylene vials feature lightweight walls for increased counting efficiency. They are premium vials, engineered and produced to precise tolerances to avoid jamming. Dimensions 27mm x 61mm (with cap on). Screw cap size 22mm. Capacity 20mL.

Catalogue No	Legacy No	Alt. No	Closure	Cap style	Inner pack qty	Pack qty
10264472	-	986701	Polypropylene	Metal foil	100	500
11556325	-	986702	Urea	Metal foil	100	500
11738839	-	986704	Polyethylene	Linerless	100	500