

Vials

accessories



Vial caps, colour coders, Fisherbrand®



- Five colours available
- Flat top allows for handwritten sample identification
- Inserts available for both internally and externally threaded vials

Cap inserts for colour coding cryogenic vials.

Catalogue No	Legacy No	Colour	For use with	Pack qty
11927974	-	White	Externally threaded vials	500
11937974	-	Yellow	Externally threaded vials	500
11947974	-	Green	Externally threaded vials	500
11957974	-	Blue	Externally threaded vials	500
11967974	-	Red	Externally threaded vials	500
11977974	-	White	Internally threaded vials	500
11987974	-	Yellow	Internally threaded vials	500
11997974	-	Green	Internally threaded vials	500
11907984	-	Blue	Internally threaded vials	500
11917984	-	Red	Internally threaded vials	500

For cryogenic vials,
please see pages 472
to 476



Vial crimpers, E-Z Crimpers™



E-Z Crimpers™ feature a cushioned ergonomic handle that is comfortable to use and can reduce hand fatigue. Crimper sizes are labelled for quick identification. Other features include polished crimping jaws for consistent sealing, a short handle and spring assembly that further aids in reducing hand fatigue. Works with the Wheaton® Crimpmaster® crimping station. An adjustable stop allows you to determine the tightness of the seal while maintaining consistent crimping.

Can be autoclaved at 121°C. Autoclaving may change the appearance and/or reduce product life.

Catalogue No	Legacy No	Alt. No	Description
10432722	-	W225303	20mm
10495951	-	W225302	13mm
10670884	-	W225301	11mm
11591215	-	W225300	8mm