

Water sampling dippers



Dippers, long handled, HDPE



For easy sampling from tanks, vats, lakes etc. Will not contaminate solutions or corrode. 29mm diameter handle is welded to the rigid bowl at a 45° angle. 3,660mm handle is available in one-piece or two-piece, screw-together construction.

Catalogue No	Legacy No	Alt. No	Capacity, mL	Handle length, mm	Handle
11302214	-	F36780-0016	500 (bowl)	910	1 piece
11312214	-	F36781-0016	500 (bowl)	1,830	1 piece
11322214	-	F36782-0016	500 (bowl)	3,660	2 piece
11332214	-	F36782-0116	500 (bowl)	3,660	1 piece
11342214	-	F36780-0032	1,000 (bowl)	910	1 piece
11352214	-	F36781-0032	1,000 (bowl)	1,830	1 piece
11362214	-	F36782-0032	1,000 (bowl)	3,660	2 piece
11372214	-	F36782-0132	1,000 (bowl)	3,660	1 piece



Drum and tank sampler, universal, Norwell



A unique product for clean and accurate sampling of drums, tanks and vats. Ideal for quality control of incoming raw materials and in-line process sampling. The product is easy to clean, with a translucent body for sample viewing and designed for ease of use. It is made from polypropylene, has a wide chemical resistance and is non-contaminating to ensure sample integrity.

An outstanding sampler for free flowing liquids. Open the safety top and the push/seal bottom, insert into the drum/tank and allow sample to flow into the tube. To close, simply push bottom against the base of the drum/tank and the seal snaps closed. The sampler can then be withdrawn.

Catalogue No	Legacy No	Alt. No	Material	Length, mm	Capacity, mL
12413060	-	D8653276	Polypropylene	1,016	300



Containers, water sampling, Thermo Scientific Sterilin™



Thermo Scientific Sterilin™ disposable water sampling containers are a convenient way of sampling water for microbiological analysis. The bottles are available pre-dosed with sodium thiosulfate or empty. All bottles have tamper-evident caps, are labelled with batch details and, where applicable, an expiry date. Bottles are sterilised by gamma radiation.

Catalogue No	Legacy No	Alt. No	Description	Capacity, mL	Pack qty
11309243	-	500WSC	Pre-dosed	500	70
11319243	-	500WSCNT	Undosed	500	70

For Fisherbrand® water sampling bottles with peel away sterilisation certificates, please see page 196.

