



General Safety

Spillage - Surface



Prevention mats, non-slip, Dychem®

04

Non-slip mats and reels have thousands of applications in the laboratory, hospital and industry. Manufactured from a specially developed polymer compound, with an exceptionally high coefficient of friction, Dychem® grips, protects, anchors and cushions. Available in a variety of types and colours, is colour-fast, odourless and resistant to most chemicals and solvents. Aids security where vibration or slipping can cause the spillage of hazardous substances or breakages. Ideal for stabilising equipment such as shakers, centrifuges, etc and for anchoring mortars, glassware, keyboards. Provides invaluable stability and shock absorption for long term experiments, as well as improving the isolation of sensitive analytical instrumentation, such as balances and microscopes, from disturbing vibrations. Available as circular or rectangular mats and in reels, may be easily cleaned with a damp cloth, or sterilised as required (do not autoclave). When wet, Dychem® is slippery, but once completely dry and clean, its full, non-slip effectiveness is restored.

Round

Catalogue No	Alt. No	Colour	Diameter, mm	Pack qty
SAT-850-010K	R211R	Red	140	2
SAT-850-020H	R211B	Blue	140	2
SAT-850-040B	R212R	Red	190	2
SAT-850-050V	R212B	Blue	190	2

Rectangular

Catalogue No	Alt. No	Colour	Width, mm	Depth, mm
SAT-860-010D	R203R	Red	350	250
SAT-860-020A	R203B	Blue	350	250
SAT-860-040R	R215R	Red	380	450
SAT-860-050X	R215B	Blue	380	450

Reels

Catalogue No	Alt. No	Colour	Width, mm	Length, mm
SAT-870-010T	R303R	Red	200	2,000
SAT-870-020Q	R303B	Blue	200	2,000
SAT-870-040K	R303W	White	200	2,000
SAT-870-050H	R304R	Red	400	2,000
SAT-870-060E	R304B	Blue	400	2,000
SAT-870-080V	R304C	White	400	2,000

2



Storage cabinets, visible, Viewcab

03

The range of Viewcab sheet steel storage cabinets are suitable for a number of working environments including laboratories, schools, universities, factories etc. The cabinets are available in three different sizes. Viewcab doors are reinforced and use 6mm laminated safety glass to BS 6206, so you can view your stored items without compromising security. Chrome lever locking handle with semi-concealed internal mechanism. Steel shelves capacity 75kg per shelf UDL. Laminated glass shelves 50kg per shelf UDL.

Mini Viewcab

Catalogue No	No. of shelves	No. of standard Winchesters held	Dimensions [w x d x h], mm
1 SAT-481-010V	2 steel shelves	18	460 x 460 x 712

Tall split Viewcab

Catalogue No	No. of shelves	No. of standard Winchesters held	Dimensions [w x d x h], mm
2 SAT-481-020S	Top unit 2 glass shelves Base unit 2 steel shelves	72	915 x 460 x 1,830

Wall mounted Viewcab

Catalogue No	No. of shelves	No. of standard Winchesters held	Dimensions [w x d x h], mm
3 SAT-481-030P	2 steel shelves	24	915 x 300 x 712

3



1

