



General Safety

Storage - Cabinets



Storage cabinets, corrosive liquid/acid

04



Manufactured from inert high density polyethylene, including hinges, these cabinets resist chemical attack better than cabinets made from any other material. Both models have facilities for padlocks and removable drip trays.

Catalogue No	Description	Dimensions [w x d x h], mm
SAT-490-031G	Central partition, two adjustable height shelves with pockets in doors for storage of small bottles	914 x 597 x 889
SAT-490-051A	Suitable for bench mounting, door can be left or right hung	362 x 413 x 495



Storage cabinets, under counter, corrosive

04



Safety features:

- Haz-Alert™ trilingual reflective labels
- U-Loc™ padlockable handle
- SpillSlope™ safety shelf
- 3 point stainless steel bullet latches

83L, 2 door, manual or self-close, blue steel safety cabinet.

Cabinet is sized and fitted with a 89mm height x 64mm depth kick plate to mount flush with cabinetry. Dual vents with flash arresters are located on the back wall for easy connection if ventilation is needed. Internal levelling feet adjust on 28.5mm centres for alignment on uneven floors. Cabinet height raises 15.9mm with leg levers fully extended. Fire resistant, double wall construction is suitable for flammable corrosives.

Not to be used with phenol, nitric acid, or sulfuric acid.

Catalogue No	Alt. No	Description	Dimensions [w x d x h], mm	No. of shelves
SAT-865-020Q	892302	Manual	889 x 559 x 889	1 (adjustable)
SAT-865-040K	892322	Self close	889 x 559 x 889	1 (adjustable)



Storage cabinet, poisons

50



Secure, robust steel cabinet for safe storage of poisonous chemicals. Constructed from epoxy resin coated galvanised steel. Can be used as a wall mounted cabinet, fittings supplied, for inside a general chemicals cabinet or placed on a bench. With a shelf. Lockable with concealed hinges. Finished in light blue/grey colour with safety label markings.

Catalogue No	Dimensions, internal [w x d x h], mm	Dimensions, external [w x d x h], mm
SAT-470-030W	300 x 287 x 327	350 x 322 x 370



Storage cabinet, poisons

03



Wall-mounted, constructed of sheet steel and finished in stoved enamel with three adjustable shelves and lockable door. The door, which has concealed hinges, is recessed for added security.

Catalogue No	Width [internal], mm	Depth [internal], mm	Height [internal], mm
SYL-400-L	450	155	605