



### Wipes, cleanroom, ISO Class 5, AlphaWipe®

04



**NEW**

ITW Texwipe's cut-edge knit wipers are used in critical environments where contamination control requirements are paramount. They are ideally suited for cleaning and wiping applications in ISO Class 5 cleanroom environments.

Constructed of 100% polyester or nylon, these wipers offer low particle release and low levels of extractables. Their superior strength allows for rigorous wiping with minimal fibre release. Premium, cut-edge, heavy-weight, all purpose polyester wiper for critical environments.

Catalogue No	Alt. No	Dimensions, mm	Pack qty
CMH-830-020R	TX1004	100 x 100 (nominal)	300
CMH-830-040L	TX1009	230 x 230 (nominal)	1,500



### Wipes, cleanroom, Class 10, ISO Class 4, Super Polx® 1200

04



- Continuous filament, non-run knitted polyester
- Laundered and packaged in a Class 10 cleanroom
- Extremely low in particle and fibre generation
- Low extractables
- Manufactured and tested with strict adherence to quality assurance guidelines
- Vacuum sealed

Catalogue No	Dimensions, mm	Pack qty
CMH-600-090K	230 x 230	150
CMH-600-100K	300 x 300	75



### Wipes, cleanroom, pre-saturated, SatPax®

04



- Resolves safety issues associated with storing, mixing, and dispensing solvents
- Semiconductor grade 2-propanol is metered for consistency
- Reverse osmosis/deionised water filtered to 0.1µm
- Resealable flexible pouches
- Control of the manufacturing process allows product customisation to meet specific cleaning criteria  
SatPax® 1000, ISO Class 5
- Cellulose/polyester

#### SatPax® 1000

Catalogue No	Solution [IPA/ H <sub>2</sub> O]	Dimensions, mm	Pack qty
CMH-600-110H	70/30	229 x 229	75
CMH-600-120E	09/91	229 x 229	75

