



## Contamination Control Hand Protection

Gloves, nitrile, powder free, BioClean Class 100,  
ISO Class 5 cleanroom compatible

04

**S**



- Hand specific cleanroom glove, processed to ensure Class 100 compatibility
- Unique packaging eliminates contamination experienced with paper use
- Pair packed in a plastic inner wrapper with a Tyvek® and polyethylene sealed outer pouch to maintain sterility
- Extended washing in pure de-ionised water ensures very low levels of particulates and extractables
- Anatomically shaped for comfort and textured for improved grip
- Latex free which eliminates allergic reactions associated with natural latex

Conforms to ISO 9002 and EN 46002. Gamma irradiated. FDA registered. Sampling and Inspection in accordance with ISO 2859. Lot numbered on inner and outer dispenser box to provide traceability and bar-coded on sealed inner box and outer carton to allow automatic data entry into stock control systems. Manufactured to exceed requirements of European Medical Glove Standard EN 455 (Part 1&2). AQL 1.5. Length 305mm. Supplied in pairs.

Catalogue No	Size	Pack type	Pack qty	Price
CMK-200-010N	5.5	Pair	50	148.79
CMK-200-020K	6	Pair	50	148.79
CMK-200-030H	7	Pair	50	148.79
CMK-200-040E	8	Pair	50	148.79
CMK-200-050B	9	Pair	50	148.79

Gloves, nitrile, powder free, surgeons,  
Profile 3000™ hand specific

04

**S**



Hand specific glove manufactured to exceed requirements of European Medical Glove standard EN 455 (1 & 2). Conforms to ISO 9002 & EN 46002. Gamma irradiated, bar-coded and lot numbered for traceability. Pairs individually packaged in a sealed pouch, for use in controlled environments. Anatomically shaped and textured for excellent grip, particularly when wet. Latex free, eliminating allergy problem.

Catalogue No	Size	Pack type	Pack qty	Price
CMK-240-020F	6	Pair	50	91.11
CMK-240-030C	6.5	Pair	50	91.11
CMK-240-040W	7	Pair	50	91.11
CMK-240-050T	7.5	Pair	50	91.11
CMK-240-060Q	8	Pair	50	91.11
CMK-240-070N	8.5	Pair	50	91.11
CMK-240-080K	9	Pair	50	91.11
CMK-240-090H	10	Pair	50	91.11

Gloves, nitrile, powder free, non-sterile,  
medical, BioClean

04



Medical grade examination glove giving extra length to ensure positive seal at wrist. Rolled cuff to prevent slipping. Suitable for clean areas, laboratories and semi controlled environments. Manufactured to exceed requirement of European glove standard EN 455 (1 & 2). Tested to AQL 1.5 level. Latex free, eliminating latex allergy problems. Ambidextrous. Length 305mm.

Catalogue No	Size	Pack qty	Price
CMK-244-040V	Small	100	36.17
CMK-244-050S	Medium	100	36.17
CMK-244-060P	Large	100	36.17
CMK-244-070M	XL	100	36.17

Gloves, nitrile, powder free, non-sterile,  
BioClean Class 100, ISO Class 5,  
cleanroom compatible, ambidextrous

**Ansell**

04



Double bagged gloves for critical environment use, thoroughly washed in deionised water for Class 100 compatibility. Latex free, eliminating latex allergy problems, textured for improved grip with an excellent tack for the safe handling of products, even when wet. Resistant to a wide range of chemicals, stronger than natural latex and anti-static. Ambidextrous. Length 305mm.

Catalogue No	Size	Pack qty	Price
CMK-260-050F	Medium	100	40.50
CMK-260-060C	Large	100	40.50
CMK-260-070W	XL	100	40.50