



Horizontal gel unit, Mini

16+

- Designed for electrophoresis with complete safety, no power supply required
- Can be run in any location and placed onto a transilluminator for visualisation during electrophoresis
- Comprises battery pack with integral lid; transparent tank; 2 x gel casting gate; 1 x 1.0mm 8 sample comb; 2 x PP3 alkaline batteries
- Gel dimensions (w x d), mm: 60 x 75
- Overall dimensions (w x d x h), mm: 110 x 170 x 30
- Maximum number of samples: 32

Catalogue No
FB57001

Accessories

Accessories	Description
FB57002	Comb 1.0mm, 4 sample
FB57004	Comb 1.0mm, 8 sample
FB57006	Comb 1.0mm, 12 sample
FB57012	Comb 1.0mm, 16 sample
FB57003	Gel casting gates, pack of 2

Gas jars, glass

- With base and ground top flange
- Without cover

Catalogue No	O.D. mm	Height, mm
GFD-205-010N	50	160
GFD-205-030H	50	210
GFD-205-050B	50	265
GFD-205-070S	80	320

Covers, glass, circular, ground one side

Catalogue No	Dia., mm	For jar dia., mm	Pack qty
GFD-220-502Y	75	50	12
GFD-220-512F	115	80	12



Heat of combustion apparatus

KS 4 16+

- Comprises borosilicate glass chamber open at the base with integral water jacket open at the top (170mm x 85mm); heat exchanger fitted to the top of the combustion chamber; wire stirrer; heat resisting platform with oxygen inlet tube and support rod; nickle crucible (15mm x 24mm) for combustion of solids; glass burner with wick and extinguishing cap for liquid fuels
- Requires retort stand for mounting
- Supplied with instructions

Catalogue No
THC-530-R

Spare Parts

Catalogue No	Description
THC-534-010N	Burner, with wick and cap
THC-534-015D	Wicks for use with THC-534-010N, pack of 5
YTL-290-030H	Nickle crucible
YTL-290-070S	Combustion chamber and water jacket

Accessories
THL-630-070N

Thermometer, -5 to 50°C



Pneumatic troughs, gas collection, glass

- Circular, for collection of gases over water
- Also suitable for storage of specimens

Catalogue No	Dia., mm	Height, mm
GFD-270-030L	200	100
GFD-270-050F	250	125
GFD-270-070W	300	125



Gas wash bottles

Drechsel (gas washing) bottle, without head

Catalogue No	Cap., mL	Socket size
QWF-350-V	125	24/29
QWF-360-X	250	24/29
QWF-370-H	500	24/29