



## Understanding Cholesterol poster

**KS 3** **KS 4** **16+**

- Dimensions (w x h), mm:  
500 x 670

Catalogue No  
ZRJ-230-060B



## Eating Disorders DVD

**NEW** **KS 4**

- This programme examines what eating disorders are, their symptoms, myths and most effective treatments
- Includes interviews with students who have suffered eating disorders, supporters of eating disorder victims, and expert psychologist Reinhold Robertson, a specialist in eating disorders
- Includes teachers' notes
- Running time: 28min

Catalogue No.  
DVD-200-180Q



## Nutrition posters

**KS 3** **KS 4** **16+**

- Set of four posters covering:  
food processing; learn to read the label; eating right; the building blocks for health



Catalogue No  
ZKH-520-020B

## Food calorimeter

**KS 4** **16+**

- Designed to determine the energy value of food and fuels
- Comprises borosilicate glass vessel (180mm x 85mm) with combustion chamber and water jacket; copper heat exchanger; clear plastic lid with wire stirrer and a hole for a thermometer (Not included); heat resistant platform with support rod, oxygen inlet tube and igniter coil mounting; 1m insulated leads
- Supplied with instructions but without stand
- Electrical supply: 6V AC/DC (Not included)
- NOTE: Oxygen supply required for complete combustion of foods

Catalogue No  
YTL-270-H

Spare Parts  
Catalogue No  
YTL-290-030H  
YTL-290-070S  
YTL-290-050B

Accessories  
Catalogue No  
PXY-400-U  
THC-534-010N  
THL-630-070N  
EKP-141-010E

Description  
Nickel crucible  
Combustion chamber and water jacket  
Igniter coils, pack of 12

Description  
Filter pump  
Burner, with wick and cap  
Thermometer, -5 to 50°C  
Multitap transformer

