

Balances

Portable balances, SG series

2 YEAR
GUARANTEE

- Sealed easy to clean keypad
- Stainless steel weighing pan
- Automatic calibration with external weight
- Capable of displaying multiple modes: g, carat, pound, ounces, troy ounces, tael - Hong Kong, tael - Singapore, tael - Taiwan, grain, pennyweight, tael - China, pound/ounces and Newton
- Applications include unit conversion, parts counting, % weighing, display hold, totalising and specific gravity function
- Supplied with power adapter
- Large easy to read high contrast LCD
- Security eyelet
- Battery saving sleep mode
- Bi-directional RS-232 option (Factory fitted) plus interface cable (Available on request)
- USB option (Factory fitted) plus interface cable and software (Available on request)



Fisherbrand

Catalogue No	Capacity, g	Readability, g	Linearity, g	Pan dimensions (dia.)/(w x d), mm
FB68975	200	0.01	±0.01	100
FB68976	400	0.01	±0.01	140 x 25
FB68977	600	0.01	±0.01	140 x 25
FB68982	600	0.1	±0.1	140 x 25
FB68985	2,000	0.1	±1	140 x 25
FB68986	5,000	1	±0.1	140 x 25
FB68987	5,000	0.1	±1	140 x 25
FB68988	10,000	1	±1	140 x 25

Accessories and Spares

Catalogue No	Description
BMT-680-090H	Alkaline battery 9V (PP3)
FB51647	Printer (including cable), with a 12 month warranty

Portable balances, Education line

NEW

- Large easy to read LCD display
- Recall function of last weighing result
- Newton Factor, percent weighing
- In use cover supplied
- Dimensions (w x d x h), mm: 194 x 225 x 67



Catalogue No	Capacity, g	Readability, g (dia.), mm	Pan size
BFM-115-010H	220	0.01	120
BFM-115-020E	220	0.01	120
BFM-115-030B	620	0.01	120
BFM-115-040V	620	0.01	120
BFM-115-050S	2,200	0.1	160
BFM-115-060P	2,200	0.1	160
BFM-115-070M	4,200	0.1	160
BFM-115-080J	4,200	0.1	160
BFM-115-090G	6,200	0.1	160
BFM-115-100G	6,200	1	160
BFM-115-110D	8,200	1	160
BFM-115-120A	8,200	1	160

† with RS232

Catalogue No	Description
Accessories	
BFE-900-514J	Auxiliary display RS/LC
BFM-110-500U	Auxiliary display with switch box
BFL-200-110T	Anti-theft device (steel cable and lock)
BFM-110-510R	Draft shield, cylindrical
BFM-110-520X	Draft shield, cylindrical, sliding doors
BFL-205-500W	In-use cover
BFM-515-530S	Lab X Direct software for simple data transfer (HB43)
BFL-205-502S	Transport case
BFM-515-540P	USB-232 converter cable
BFM-515-545F	Weigh pan, 120mm