



Compression spring balance

- Capacity: 500g x 5g
- Dimensions, mm: 175 x 145 x 40

Catalogue No
BAW-560-U



Weighing funnels, borosilicate glass

Fisherbrand®

- Ideal for weighing small quantities of powder
- Powder can be introduced via the tubular stem

Catalogue No	Capacity, mL	Length, mm
FB51970	3	70
FB51971	6	85
FB51972	10	100



Weighing bottles, glass

- Supplied in packs of 4

Catalogue No	Capacity, mL	Height, mm	Diameter (external), mm
Tall, stoppered			
BGW-600-010Y	8	40	20
BGW-600-030C	15	50	25
BGW-600-050T	28	60	30
BGW-600-070N	60	80	40
BGW-600-090H	135	100	50
Tall, capped			
BGW-610-010B	8	40	20
BGW-610-030S	15	50	25
BGW-610-050M	28	60	30
Squat, stoppered			
BGW-620-010R	30	25	50
BGW-620-030L	55	30	60

Platform type balance

- With zero adjustment
- Dimensions (w x d x h), mm: 285 x 275 x 45
- Capacity: 125g x 1kg

Catalogue No
BAW-810-Y



Platform force meter, Physics,
Measurement, page 7
Spring balances, Physics,
Measurement, page 6



Weighing boats, disposable, polystyrene

Fisherbrand®

- Anti-static range alleviates powder hang up
- Flexible to enable ease of transfer of contents
- Flat bases for safe and firm positioning on balance pan
- Smooth radiused corners and sloping side for complete sample transfer
- Economical

Catalogue No	Style	Cap., mL	Dimensions (w x d), mm	Pack qty
Standard white				
FB50323	Diamond	5	55 x 30	1,000
FB50325	Diamond	30	80 x 60	500
FB50327	Diamond	100	120 x 100	250
FB50333	Square	7	44 x 44	1,000
FB50335	Square	100	80 x 80	1,000
FB50337	Square	250	140 x 140	1,000
Anti-static white				
FB50343	Diamond	5	55 x 30	1,000
FB50345	Diamond	30	80 x 60	500
FB50347	Diamond	100	120 x 100	250
FB50357	Square	250	140 x 140	1,000



Weighing bottles, polypropylene



- Lightweight, reusable translucent bottles with airtight push on lids
- Highly resistant to acid attack
- Conforms to food and agricultural regulations: EEC 90/128 92/93 and FDA 21CFR 177.1520
- Supplied in packs of 10

Catalogue No	Cap., mL	Height, mm	Tare mass, g
BGW-320-010S	25	30	8
BGW-320-030M	30	50	8
BGW-320-050G	50	30	10
BGW-320-070A	65	60	14
BGW-320-090R	65	35	18
BGW-320-110X	200	90	28