

Evaporating basins, porcelain, flat bottom with spout

- Shallow form
- Glazed except outside base

Catalogue No	Cap., mL	Dia., mm	Height, mm
BHW-300-010V	20	63	13
BHW-300-030P	45	80	16
BHW-300-070D	100	100	20
BHW-300-210N	1,100	200	55



Evaporating basins, porcelain, round bottom with spout

- Shallow form
- Sizes up to 285mL glazed except rim, 420mL and larger glazed except outside base

Catalogue No	Cap., mL	Dia., mm	Height, mm
BHW-380-030F	50	75	30
BHW-380-050W	75	80	32
BHW-380-070Q	95	85	36
BHW-380-090K	115	100	40
BHW-380-130B	225	110	45
BHW-380-150S	285	125	50
BHW-380-170M	420	150	45
BHW-380-190G	580	160	64
BHW-380-210D	1,000	200	80
BHW-380-250X	2,000	250	100
BHW-380-270Y	3,000	315	126



Evaporating basins, porcelain, flat bottom with spout

- Deep form, French shape
- Glazed inside and outside

Catalogue No	Cap., mL	Dia., mm	Height, mm
BHW-320-030B	30	58	24
BHW-320-050S	60	75	30
BHW-320-070M	100	84	38
BHW-320-090G	140	97	40
BHW-320-110D	250	110	50
BHW-320-130U	350	125	52
BHW-320-150X	500	150	60



Evaporating basins, porcelain, round bottom with spout

- Deep form, French shape
- Glazed inside and out

Catalogue No	Cap., mL	Dia., mm	Height, mm
BHW-420-030Q	50	75	40
BHW-420-110S	350	125	65



Incinerating basins, porcelain, rectangular with conical walls

- Glazed except outside base

Catalogue No	Width, mm	Depth, mm	Height, mm
BHW-550-030K	48	25	10
BHW-550-050E	55	40	15
BHW-550-070V	70	45	12
BHW-550-110M	95	58	15

