

Bottles



Reagent bottles, glass, narrow neck

- Interchangeable octagonal head, polypropylene stopper
- Available in clear or amber glass

Catalogue No	Cap., mL	Neck type	Height, mm	O.D., mm	Pack qty
Clear					
BTF-201-030J	30	16/16	85	35	10
BTF-201-040G	50	16/16	100	42	10
BTF-201-071S	100	16/16	120	57	10
BTF-201-090X	250	24/20	145	65	10
BTF-201-111J	500	24/20	177	82	2
BTF-201-131D	1,000	29/32	226	102	2

Amber					
Catalogue No	Cap., mL	Neck type	Height, mm	O.D., mm	Pack qty
BTF-221-051K	50	16/16	100	42	10
BTF-221-071E	100	16/16	120	57	10
BTF-221-091V	250	24/20	145	65	10
BTF-221-110U	500	24/20	177	82	2
BTF-221-130X	1,000	29/32	226	102	2



Reagent bottles, PVC coated glass, narrow neck

- 250mL capacity bottles with labels and PVC coating
- Labels have hazard symbols and phrases where required
- Interchangeable polypropylene hexagonal head stoppers
- Can be dry heat sterilised at 140°C, ensure PVC surfaces are not touching
- Prevent immersion for prolonged periods in cleaning fluids, as PVC coating may be attacked
- Do not autoclave
- Supplied in packs of 2

Catalogue No	Label	Stopper
Clear		
BTJ-550-030Y	Acetic acid, dil.	Plastic
BTJ-550-070T	Ammonia solution, dil.	Plastic
BTJ-550-090N	Ammonia, conc.	Plastic
BTJ-550-170P	Ethanol	Plastic
BTJ-550-230A	Hydrochloric acid, dil.	Plastic
BTJ-550-250R	Hydrochloric acid, conc.	Glass
BTJ-550-290F	Lime water	Plastic
BTJ-550-310C	Nitric acid, dil.	Plastic
BTJ-550-430P	Propanone (acetone)	Plastic
Amber		
BTJ-550-130E	Bromine water	Glass
BTJ-550-330T	Nitric acid, conc.	Glass
BTJ-550-450J	Silver nitrate	Plastic



Fisherbrand

Reagent and media bottles, borosilicate clear glass, wide screw neck

- Chemically resistant graduated bottles supplied with blue polypropylene cap and pouring ring
- ISO/GL 45 screw thread standard to all sizes
- Conforms to ISO R115 and DIN 168
- Can be repeatedly sterilised to 140°C
- Do not autoclave bottles with a tightly screwed cap

Catalogue No	Cap., mL	Height, mm	O.D., mm	Pack qty
FB33144	100	105	56	10
FB33145	250	143	70	10
FB33146	500	181	86	10
FB33147	1,000	230	101	10
FB33148	2,000	267	136	each

Polypropylene screw caps, with lip seal, autoclavable to 140°C, packs of 10

Catalogue No	ISO thread GL	Colour
FB50487	45	Blue
FB50488	45	Green
FB50489	45	Yellow
FB50490	45	Black

Polypropylene pouring rings, autoclavable to 140°C, packs of 10

FB50491	45	White
FB50492	45	Red
FB50493	45	Blue
FB50494	45	Green
FB50495	45	Yellow
FB50496	45	Black



Winchester bottles, glass, narrow screw neck

- Supplied with black plastic screw cap and woodpulp faced steran liner (PVDC)
- Available in clear or amber glass

Catalogue No	Cap., mL	Thread style	Pack qty
Clear			
BTF-430-010H	30	R3/20	40
BTF-430-055Y	100	R3/28	20
BTF-430-080J	250	R3/28	20
BTF-430-110D	500	R6/31	10
BTF-430-130U	1,000	R6/31	10
Amber			
BTF-440-005Q	15	R3/20	40
BTF-440-010A	30	R3/20	40
BTF-440-040X	50	R3/24	40
BTF-440-050L	100	R3/24	20
BTF-440-045E	120	R3/24	20
BTF-440-070F	150	R3/24	20
BTF-440-080C	200	R3/24	20
BTF-440-110T	500	R6/31	10
BTF-440-130N	1,000	R6/31	10
BTF-440-150H	2,000	R6/31	10
BTF-440-140K	2,500	R6/31	12