



Aspirator bottles, polyethylene, narrow screw neck

- 8mm screw in tap and screw on cap with plastic covered, hinged handle

Catalogue No	Cap., mL
ASP-520-010V	5,000
ASP-520-030P	10,000

Spare Part Catalogue No	Description
ASP-524-010U	Screw in tap



Aspirator bottles, rectangular, polyethylene, narrow screw neck

- 15mm tap, screw on cap and carrying handle

Catalogue No	Capacity, mL
ASP-560-010Q	25,000

Spare Part Catalogue No	Description
ASP-564-010P	Tap



Aspirator bottles, heavy duty HDPE, wide screw neck

- Sturdy construction having a 2mm wall thickness whilst retaining translucency
- $\frac{3}{4}$ " BSP Waddington Duval tap mates with a 12mm deep aspirator for total leakproof performance
- Recommended that the tap is greased before assembly

Catalogue No	Cap., mL	Neck size, mm	Height, mm	O.D. mm
ASP-650-010P	5,000	45	335	165
ASP-650-030J	10,000	53	420	210
ASP-650-050D	25,000	53	540	275
ASP-650-070U	50,000	53	660	365

Spare Part Catalogue No	Description
ASP-654-010X	Tap



Aspirator bottle, rectangular, HDPE, wide screw neck

- Designed to lie on it's side so that tap overhangs the edge of a table or shelf
- Rounded corners, sturdy handle and bottom grip for pouring
- $\frac{1}{4}$ turning type tap which extends to 108mm, with a 35mm long outlet tube, 16mm in diameter
- 52mm diameter fill opening with screw cap and 26mm side opening with a screw cap


Catalogue No	Cap., mL	Height, mm	Width, mm	Depth, mm
ASP-890-010R	20,000	380	250	250



General purpose bottles, polystyrene, wide neck

- White urea caps with steran faced woodpulp wad

Catalogue No	Cap., mL	Neck size, mm	Height, mm	O.D., mm	Pack qty
BTK-352-030C	30	26	48	35	20
BTK-352-050T	60	32	59	43	20
BTK-352-070N	125	36	68	54	20
BTK-352-090H	250	54	79	69	10
BTK-352-110E	350	64	87	80	10



Aspirator bottle rack

- Manufactured from, corrosion-resistant plastic material shipped in a compact package ready for assembly
- The five-shelf rack has 16 leg components that slide into the shelf corner sockets

- Legs are 38mm in diameter, shelf spacing is 260mm
- Dimensions between legs are 270mm x 310mm
- Each rack can support 91kg or 18kg per shelf
- Shipping weight is 4.7kg

Catalogue No
ASP-890-500A



General purpose bottles, medium weight polyethylene, extra wide neck

- Self sealing polypropylene cap

Catalogue No	Cap., mL	Height, mm	O.D., mm	Pack qty
BTK-420-070P	150	68	60	10
BTK-420-090J	250	82	78	10
BTK-420-110G	500	86	117	5
BTK-420-130A	1,000	148	117	5