

Bottles / Brushes

Wash bottles, round LDPE, narrow screw neck

- Sloping shoulder and interchangeable leakproof closures
- Printed with a choice of water or distilled water
- Supplied in packs of 5
- Bottles marked acetone, methanol and other solvents have been withdrawn on the advice of CLEAPSS

These are available as specials, please contact 01509 233344 for details

Catalogue No	Capacity, mL	Wording	Print and closure
BTS-640-030J	250	Water	White
BTS-640-150W	250	Distilled water	White
BTS-640-210H	500	Water	White
BTS-640-330U	500	Distilled water	White



Wash bottle carrier, polypropylene



- Capable of carrying up to six wash bottles
- Holds either 250mL or 500mL bottles
- Good chemical resistance
- Dimensions (w x d x h), mm: 286 x 194 x 225

Catalogue No
BTS-920-030W



Utility bottle carrier, polypropylene



- Accommodates square bottles and rectangular bottles and unwire test tube racks
- External dimensions (w x d x h), mm: 375 x 241 x 114
- Internal dimensions (w x d x h), mm: 89 x 340 x 111

Catalogue No
BTS-844-030P



Balance dusting brush, hair

- Double ended of pure blended hair
- Suitable for dusting laboratory balances, cleaning other instruments and weighing scoops
- Large brush sweeps large areas, small brush cleans crevices

Catalogue No	Visible length of hair, mm	Overall length, mm
BUR-210-B	Large = 25, Small = 10	145



Wash bottles, polyethylene, wide screw neck

- Designed to minimise rippling and give an accurate controlled stream of liquid
- Colour coded closures for easy identification
- Graduated
- Supplied in packs of 5

Catalogue No	Capacity, mL	Neck size, mm	O.D., mm	Closure/stem/ wording colour
BTS-500-030W	250	38	60	Neutral
BTS-500-050Q	250	38	60	Blue
BTS-500-070K	250	38	60	Yellow
BTS-500-090E	250	38	60	Red
BTS-500-110B	500	38	75	Neutral
BTS-500-130S	500	38	75	Blue
BTS-500-150M	500	38	75	Yellow
BTS-500-170G	500	38	75	Red



Winchester bottle carrier

- Suitable for carrying plain Winchester type bottles 2.5L capacity

Catalogue No	Type
BTS-865-030S	Single
BTS-865-050M	Double



Bottle carrier, 2 place

- Polyethylene coated mild steel carrier for two 2.5L bottles of either Winchester pattern or Eurobottle type

- Maximum diameter of bottle held in 150mm

Catalogue No
BTS-870-E

Winchester bottle carrier

- Accommodates one or two 2.5L Winchester shape bottles
- Liquid tight, with carrying handle pivoted for easy pouring

Catalogue No	Type	Colour
BTS-850-S	Single	Red
BTS-852-G	Double	White



Pencil brush, hair

- In metal mount on wood handle
- Supplied in packs of 10

Catalogue No	Head dia., mm	Length of head, mm
BUR-860-K	8	18