


Burettes



Burettes, glass, single bore, pinchvalve, BS846 Class B

- Ideal burettes for educational use
- Graduated
- Detachable rubber pinch valve and separate glass jet
- Pinch valves consist of rubber moulding with a bore at each end tapering from 7 to 4mm

Catalogue No	Cap., mL	Grads., mL
BWF-301-A	50	0.10

Spare Parts

Catalogue No	Description
BWF-330-Y	Burette pinchvalves, pack of 5
BWF-350-R	Burette jets, glass, pack of 5



Burettes, glass, single bore, interchangeable PTFE stopcock, BS846 Class B

- Blue ceramic graduations

Catalogue No	Cap., mL	Grads., mL	Tolerance, mL
FB50041	1	0.01	±0.02
FB50045	10	0.05	±0.05
FB50049	25	0.10	±0.10
FB50051	50	0.10	±0.10
FB50053	100	0.20	±0.20

Burettes, glass, single bore, interchangeable PTFE stopcock, BS846

Catalogue No	Cap., mL	Grads., mL
Class B		
BWF-740-050A	10	0.02
BWF-740-090L	25	0.05
BWF-740-110Y	25	0.10
BWF-740-150T	50	0.10
BWF-740-170N	100	0.20

Class B Schelbach

BWF-750-110B	25	0.10
BWF-750-150M	50	0.10
BWF-750-170G	100	0.20

Class A, uncertificated, individually tested for accuracy

BWF-780-070M	10	0.02
BWF-780-100G	25	0.05
BWF-780-130U	25	0.10
BWF-780-150X	50	0.10
BWF-780-170Y	100	0.20

Class A, with certificate of calibration

BWF-785-070F	10	0.02
BWF-785-100W	25	0.05
BWF-785-130N	25	0.10
BWF-785-150H	50	0.10
BWF-785-170B	100	0.20



Burettes, glass, single bore, interchangeable glass stopcock, BS846 Class B

Fisherbrand®

- Blue ceramic graduations

Catalogue No	Cap., mL	Grads., mL	Tolerance, mL
FB50009	25	0.10	±0.10
FB50011	50	0.10	±0.10
FB50013	100	0.20	±0.20

Burettes, borosilicate glass, single bore, interchangeable PTFE stopcock, Class B, BS-ISO385:2005

- White graduations



Catalogue No	Cap., mL	Grads., mL
BWF-470-010P	10	0.02
BWF-470-015F	25	0.10
BWF-470-020M	50	0.10
BWF-470-025C	100	0.20

