



Centrifuge, Labofuge® 200

- Compact benchtop centrifuge designed for the safe centrifugation of blood samples making it ideal for use in medical practises, small laboratories and as a stand by unit in large laboratories
- Microprocessor control ensures ease of use and reproducibility of results
- Supplied with 12 place autoclavable rotor for 15mL DIN tubes
- Maintenance free brushless induction drive

Technical Specification

| | |
|----------------------------------|-------------------|
| Max. speed, rpm | 5,300 |
| Min. speed, rpm | 1,600 |
| Max. RCF, xg | 3,030 |
| Dimensions (w x d x h), mm | 284 x 375 x 240 |
| Mass, kg | 10.7 |
| Electrical supply | 220/230V, 50/60Hz |
| Fixed angle rotor - CFH-135-010Y | |
| Max. capacity, mL | 12 x 15 |
| Max. speed, rpm | 5,300 |
| Max. RCF, xg | 3,030 |
| Acceleration time, s | 40 |
| Braking time, s | 45 |
| Max. radius, mm | 96.5 |

Catalogue No
CFH-130-010P

Rotors

| Catalogue No | Description |
|--------------|--|
| CFH-135-010Y | Fixed angle rotor, 12 x 15mL, approved for 10,000 cycles |

Accessories

| | |
|--------------|--|
| CFH-135-500X | Adapters, for 5/7mL, complete set for one rotor |
| CFH-135-502K | Pads, for 10mL Sarstedt Monovettes, complete set for one rotor |
| CFH-135-504G | Protective sleeves, complete set for one rotor |

Centrifuge tubes

| Catalogue No | Cap., mL | Material | Grads., mm | O.D., (max.) | Height, mm | RCF xg | Pack qty |
|--------------|----------|---------------|------------|--------------|------------|--------|----------|
| CFT-480-A | 15 | Soda glass | - | 17 | 111 | 1,300 | 10 |
| CFS-110-X | 12 | Pyrex | - | 17 | 111 | 3,000 | each |
| CFT-440-F | 15 | Polycarbonate | Yes | 17 | 119 | 6,000 | 10 |
| CFS-106-Y | 12 | Polycarbonate | - | 16 | 100 | 50,000 | 10 |
| CFT-450-V | 15 | Polycarbonate | - | 16 | 116 | 50,000 | 10 |



For latest offers and information...

Why not sign up now for our e-newsletters at:

www.griffineducation.co.uk/newsletter