



Paper chromatography kit

- Comprehensive kit suitable for introducing the theory and principles of ascending, descending and reverse phase chromatography
- Comprises:
 - 1 x cylindrical tank, glass 75mm x 305mm, with ground glass cover
 - 1 x rectangular tank, glass, 230mm x 90mm x 230mm, with ground glass cover
 - 1 x semi-circular trough, glass, 205mm x 24mm
 - 2 x wire frames, nylon coated, for use in rectangular and cylindrical tanks
 - 1 x reel of chromatography paper strip, 100m x 50mm
 - 1 x box of chromatography paper sheets, 100 sheets, 200mm x 200mm
 - 1 x manual of experiments by D. R. Browning, B.Sc., F.R.I.C., A.R.T.C.S.

Catalogue No
CJA-200-L

Spare Parts

Catalogue No	Description
GFD-220-512F	Covers for cylindrical tank, pack of 12
CJA-221-010L	Glass tank, with cover
CJA-220-070W	Trough, semicircular
CJF-230-050R	Chromatography paper strips, each
CJF-250-030J	Chromatography paper sheets, pack of 100
CJA-200-805P	Rectangular wire frame
GFD-205-070S	Tank, cylindrical



Thin layer chromatography kit

- Comprehensive kit suitable for introducing the principles of adsorption chromatography using the thin layer technique
- Comprises:
 - 1 x applicator, for plate material, to accept glass plates (150mm x 75mm) and microscope slides (75mm x 25mm)
 - 5 x glass plates (150mm x 75mm)
 - 100 x microscope slides (76mm x 26mm)
 - 1 x rectangular glass tank with cover (230mm x 90mm x 230mm)
 - 1 x Coplin jar
 - 100 x capillary tubes
 - 1 x manual of experiments by D. R. Browning, B.Sc., F.R.I.C., A.R.T.C.S.

Catalogue No
CJK-300-S

Spare Parts

Catalogue No	Description
CJA-221-010L	Glass tank, with cover
YSH-102-L	Microscope slides, pack of 50
MNK-820-H	Coplin jar
MFB-210-538L	Capillary tubes, open both ends, 75mm, pack of 100
Accessories	
FPM-620-090P	Separating funnel, 250mL



Update your contact
details online...

www.griffineducation.co.uk/mydetails