

Electrical power supplies / Filtration



Power supply,
battery replacement

- Low cost power supply housed in a robust ABS plastic case suitable for replacing batteries in various applications
- Smooth and regulated DC
- Voltage outputs selected via rotary selector switch (1.5V, 3V, 4.5V or 6V at currents up to 1A)
- Two 4mm colour coded 4mm output sockets
- Dimensions (w x d x h), mm: 100 x 150 x 80
- Mass, g: 800
- Electrical supply: 230V, 50Hz/60Hz

Catalogue No
EKR-451-010W



Filter papers,
schools grade

- Supplied in packs of 100

Catalogue No	Diameter, mm
FDH-190-109J	90
FDH-190-129D	110
FDH-190-149U	125
FDH-190-169X	150
FDH-190-189Y	180



Filtration, quartz wool

A mass of fibres, average diameter 6µm, used for filtering acids and hot gases or as a plug for combustion tubes.

Catalogue No	Quantity, g
FDK-960-E	10



Power supply, mains
adapters,
BS415:1990

- Mains adapters suitable for a wide variety of low voltage equipment
- Voltage output can be selected to suit exact requirements
- LED power on indicator
- Negative polarity facility
- Supplied with multi-way plug and 2m output lead

Catalogue No	Current (max.), mA	DC volts
EKR-460-010V	300	3, 4, 5, 6, 7.5, 9, 12
EKR-460-020S	1,200	1.5, 3, 4, 5, 6, 7.5, 9, 12



Filter papers, cellulose acetate

- Provide excellent thermal stability
- Low absorption and high wet strength
- Can be in-lined steamed, autoclaved or dry heat sterilised
- Please enquire for other pore sizes

Catalogue No	Pore size, µm	Dia., mm	Pack qty
FDE-100-010G	0.2	25	100
FDE-100-020D	0.2	47	100
FDE-100-030A†	0.2	47	100
FDE-100-040U	0.2	142	25
FDE-100-050R	0.2	142	100
FDE-105-010W	0.45	25	100
FDE-105-020T	0.45	47	100
FDE-105-025T	0.45	85	100
FDE-105-030Q†	0.45	47	100
FDE-105-040N	0.45	142	25
FDE-105-050K	0.45	142	100

†Sterile