



Filter papers, qualitative, QL range **Fisherbrand**

- Suitable for general filtration of many types of samples
- Supplied in packs of 100

Catalogue No Diameter, mm

Fisherbrand QL100

May be used in many applications where Whatman No. 1 or Schleicher & Schuell No. 595 are used

Medium retention and medium/fast flow rate

FB59011	42.5
FB59015	55
FB59017	70
FB59019	85
FB59021	90
FB59023	110
FB59025	125
FB59027	150
FB59029	185
FB59031	240
FB59033	270
FB59035	320
FB59037	460 x 570 sheets

Fisherbrand QL110

May be used in many applications where Whatman No. 1 or Schleicher & Schuell No. 595 are used

Medium/fine retention and medium flow rate

FB59041	70
FB59043	125
FB59045	150
FB59047	185
FB59049	240
FB59051	270
FB59055	460 x 570 sheets

Fisherbrand QL120

May be used in many applications where Whatman No. 3 or Schleicher & Schuell No. 598 are used

Medium/fine retention and high loading capacity

FB59075	70
FB59077	90
FB59081	125
FB59083	150



Filter papers, qualitative

Catalogue No Diameter, mm Pack qty

Filter papers No. 1

- Medium retention and - medium/fast flow rate
- Widely used for air pollution analysis
- Pore size 11µm

FDH-240-035Q	30	400
FDH-240-050U	42.5	100
FDH-240-070X	55	100
FDH-240-090Y	70	100
FDH-240-100Y	85	100
FDH-240-110F	90	100
FDH-240-130W	110	100
FDH-240-150Q	125	100
FDH-240-170K	150	100
FDH-240-190E	185	100
FDH-240-210B	240	100
FDH-240-230S	270	100
FDH-240-250M	320	100
FDH-240-330X	460 x 570 sheets	100
FDH-240-350Y	580 x 680 sheets	100

Filter papers No. 2

- Slightly thicker and more retentive than No. 1
- Pore size 8µm

FDH-260-070A	55	100
FDH-260-090R	70	100
FDH-260-110X	90	100
FDH-260-130Y	110	100
FDH-260-170T	150	100
FDH-260-190N	185	100
FDH-260-210K	240	100

Filter papers No. 3

- A thicker paper, offering improved wet strength and particle retention compared with No. 1
- Suitable for use in Buchner funnels
- High precipitate capacity without clogging
- Pore size 6µm

FDH-280-070J	55	100
FDH-280-090D	70	100
FDH-280-110A	90	100
FDH-280-130R	110	100
FDH-280-150L	125	100
FDH-280-170F	150	100
FDH-280-190W	185	100
FDH-280-210T	240	100
FDH-280-330J	460 x 570 sheets	100