



Flame testing

Catalogue No Description

FGL-340-A	Blowpipe, bent, plated brass (200mm), pack of 10
FGL-520-V	Glasses, for flame test, (75mm x 50mm), cobalt blue, for the detection of potassium in the presence of sodium, etc, pack of 100
FGL-550-A	Rods, silica, for flame test, plain, (75mm x 1mm), pack of 100
PRF-440-S	Wire, platinum, for flame test, (38mm x 0.38mm), fused into glass handle (75mm)



Cork rings, flasks,

- Dimensions quoted as approximate

Catalogue No	Dia. (ext.), mm	Dia. (int.), mm	Height, mm
CSD-340-071H	195	160	45
CSD-340-091B	240	200	45
CSD-340-051N	150	120	45
CSD-340-011C	75	45	45
CSD-340-031T	115	85	45



Stabilising lead weights, flasks, vinyl coated

- Flexible lead weights coated with chemical resistant red Vikem® vinyl
- Can be coiled into a helix, spiral or cone to stabilise vessels in a water bath
- Supplied in packs of 5

Catalogue No	Length, mm	Dia. (max. ext.), mm
FHM-330-010G	300	6.4
FHM-330-030A	300	8



Flasks, Pyrex® glass, round bottom, narrow neck, BS2734, ISO1773

- With beaded rim

Catalogue No	Cap., mL	Neck I.D., mm	Height, mm	Pack qty
FHB-275-030X	100	18	115	5
FHB-275-050Y	250	30	145	5
FHB-275-070C	500	30	175	2
FHB-275-090T	1,000	36	210	2
FHB-275-110Q	2,000	45	260	each
FHB-275-120N	10,000	57	420	each