

Incubators, general purpose



The GENLAB General Purpose incubators offer a range of highly efficient, reliable, cost effective units to suit most drying, sterilising and general laboratory applications

- Natural convection on all models
- Exterior constructed from sheet steel finished in an easy clean powder coated paint
- Interior chamber made from mild steel coated with aluminium
- Thermostat control with independent safety overheat facility both with calibrated scales and tamper proof locks
- Horizontal version available for the 75L size
- Also available with digital control and viewing window options
- Fans available on 30L+ models, contact us for details
- Temperature range:
 - Ambient +5 to 100°C (Without fan)
 - Ambient +5 to 60°C with viewing window
- Fluctuation: $\pm 0.25^{\circ}\text{C}$ at 37°C



Catalogue No	Description	Cap., L	Dimensions, int. (w x d x h), mm	Dimensions, ext. (w x d x h), mm	Shelves, no. of	Shelf positions, no. of	Weight, kg
INC-400-006Y	Vertical	6	230 x 190 x 150	320 x 330 x 390	1	1	7
INC-410-065Y	Vertical, window	6	230 x 190 x 150	320 x 330 x 390	1	1	7
INC-400-018B	Vertical	18	260 x 270 x 260	380 x 410 x 520	2	2	14
INC-405-018R	Vertical, digital control	18	260 x 270 x 260	380 x 410 x 520	2	2	14
INC-410-018R	Vertical, window	18	260 x 270 x 260	380 x 410 x 520	2	2	14
INC-400-030L	Vertical	30	360 x 350 x 240	480 x 490 x 500	2	2	20
INC-405-035R	Vertical, digital control	30	360 x 350 x 240	480 x 490 x 500	2	2	20
INC-410-030E	Vertical, window	30	360 x 350 x 240	480 x 490 x 500	2	2	20
INC-400-040Y	Vertical	40	360 x 350 x 320	450 x 490 x 580	2	3	23
INC-405-045X	Vertical, digital control	40	360 x 350 x 320	450 x 490 x 580	2	3	23
INC-410-040B	Vertical, window	40	360 x 350 x 320	450 x 490 x 580	2	3	23
INC-400-050F	Vertical	50	490 x 330 x 330	610 x 470 x 590	2	3	26
INC-405-055L	Vertical, digital control	50	490 x 330 x 330	610 x 470 x 590	2	3	26
INC-410-050V	Vertical, window	50	490 x 330 x 330	610 x 470 x 590	2	3	26
INC-400-075M	Vertical	75	490 x 450 x 330	610 x 600 x 590	2	3	34
INC-405-755G	Vertical, digital control	75	490 x 450 x 330	610 x 600 x 590	2	3	34
INC-410-750Q	Vertical, window	75	490 x 450 x 330	610 x 600 x 590	2	3	34
INC-415-075V	Horizontal	75	530 x 350 x 420	850 x 500 x 560	2	3	39
INC-420-050X	Horizontal	50	360 x 350 x 420	680 x 500 x 560	2	3	31
INC-425-060E	Horizontal, digital control	50	360 x 350 x 420	680 x 500 x 560	2	3	31
INC-430-065R	Horizontal, window	50	360 x 350 x 420	680 x 500 x 560	2	3	31
INC-455-075Q	Horizontal, digital control	75	530 x 350 x 420	850 x 500 x 560	2	3	39
INC-450-075A	Horizontal, window	75	530 x 350 x 420	850 x 500 x 560	2	3	39