



Microscope, Pyser XES-50-P2 Stereo

- Simple to use
- Excellent optical quality
- 20x total magnification
- Upright binocular head with adjustable interpupillary distance from 55 to 75mm
- Paired WF 10x eyepieces with ± 5 diopter adjustment on one eyepiece
- Slide in 2x objective
- 110mm x 120mm plain stand with removable white/black plate, 94.5mm diameter and spring clips
- Rack and pinion focusing

Catalogue No
MIT-220-030L



Microscope, Pyser XES-70C-2L Stereo

- 20x and 40x total magnification
- Inclined binocular head with adjustable interpupillary distance from 55 to 75mm
- Paired WF 10x eyepieces with ± 5 diopter adjustment on one eyepiece
- Rotatable turret with choice of 2x and 4x objectives
- Illuminated stand with removable and interchangeable white/black and frosted glassplate, 94.5mm diameter with spring clips
- Incident and transmitted light sources, separately switchable
- Rack and pinion focusing

Catalogue No
MIT-230-020H



Microscope,
pocket, portable

Fisherbrand

NEW

- Ideal for inspection; field use; metal analysis; cloth/paper examination and plant studies
- Adjustable light source, adjustable focus control knob
- Chemical and shock resistant ABS plastic
- FB70286 supplied with vinyl carrying case and FB70287 supplied with a snap-on slide holder
- Dimensions (w x h), mm:
140 x 38 x 22
- Requires AA batteries (Included)

Catalogue No
FB70286
FB70287

Description
Illuminated microscope 30x with carrying case
Illuminated microscope, zoom