

## General purpose ovens

- GENLAB general purpose ovens are efficient, cost effective units to suit most drying, sterilising and general laboratory applications
- Exterior constructed from sheet steel finished in an easy clean powder coated paint
- Interior chamber is made from mild steel coated with aluminium
- Thermostat control with independent safety overheat facility both with calibrated scales and tamper proof locks
- Horizontal version available for the 50L and 75L sizes
- 30L to 75L also available with digital control
- Temperature range, °C: +40 to +450
- Fluctuation:  $\pm 0.75^{\circ}\text{C}$  @  $100^{\circ}\text{C}$



 **Dual purpose ovens/incubators, Incubators, General Labware, page 433**

| Catalogue No | Description | Cap., L | Dimensions, (w x d x h), mm |                 | Number of shelves | Number of shelf positions | Weight, kg |
|--------------|-------------|---------|-----------------------------|-----------------|-------------------|---------------------------|------------|
|              |             |         | Internal                    | External        |                   |                           |            |
| OVB-100-006U | Vertical    | 6       | 230 x 190 x 150             | 320 x 330 x 390 | 1                 | 1                         | 7          |
| OVB-100-018N | Vertical    | 18      | 260 x 270 x 260             | 380 x 410 x 520 | 2                 | 2                         | 14         |
| OVB-100-030A | Vertical    | 30      | 360 x 350 x 240             | 480 x 490 x 500 | 2                 | 2                         | 20         |
| OVB-105-035G | Vertical†   | 30      | 360 x 350 x 240             | 480 x 620 x 500 | 2                 | 2                         | 20         |
| OVB-100-040U | Vertical    | 40      | 360 x 350 x 320             | 450 x 490 x 580 | 2                 | 3                         | 23         |
| OVB-105-045D | Vertical†   | 40      | 360 x 350 x 320             | 450 x 620 x 580 | 2                 | 3                         | 23         |
| OVB-100-050R | Vertical    | 50      | 490 x 330 x 30              | 610 x 470 x 590 | 2                 | 3                         | 26         |
| OVB-105-055A | Vertical†   | 50      | 490 x 330 x 330             | 610 x 590 x 590 | 2                 | 3                         | 26         |
| OVB-100-075B | Vertical    | 75      | 490 x 450 x 330             | 610 x 600 x 590 | 2                 | 3                         | 34         |
| OVB-105-075F | Vertical†   | 75      | 490 x 450 x 330             | 610 x 720 x 590 | 2                 | 3                         | 34         |
| OVB-150-010R | Horizontal  | 50      | 360 x 350 x 420             | 680 x 500 x 560 | 2                 | 3                         | 31         |
| OVB-150-020Y | Horizontal* | 50      | 320 x 350 x 450             | 680 x 500 x 560 | 2                 | 3                         | 31         |
| OVB-150-030L | Horizontal† | 50      | 320 x 350 x 450             | 680 x 500 x 560 | 2                 | 3                         | 31         |
| OVB-150-075M | Horizontal  | 75      | 490 x 450 x 330             | 530 x 350 x 420 | 2                 | 3                         | 39         |
| OVB-155-750Q | Horizontal† | 75      | 490 x 450 x 330             | 490 x 350 x 450 | 2                 | 3                         | 39         |

† With fan and digital control

\* With fan only



## Paper towels

- General purpose towels available in embossed 1-ply or 3-ply with extra softness
- Premium quality for maximum comfort
- 100% recycled fibre, no dye, no bleach
- Made from AIRFLEX® highly absorbent fabric
- Excellent value in everyday use



## Paper wipes, blue, WYPALL® L30

- 2-ply centre feed roll
- Patented AIRFLEX® quilting gives increased absorbency and wet strength
- 125mm roll diameter
- Sheet size 240mm x 460mm
- 82 sheets per roll
- 24 rolls per pack

Catalogue No  
CMC-717-022U

| Catalogue No | Description   | Dimensions, mm | Pack qty         |
|--------------|---------------|----------------|------------------|
| CMC-710-021S | 1-ply white   | 240 x 310      | 24 sleeves x 150 |
| CMC-710-030R | 2-ply natural | 420 x 220      | 30 sleeves x 84  |

| Accessory    |                            |
|--------------|----------------------------|
| Catalogue No | Description                |
| CMC-710-500G | Dispenser for CMC-710-030R |

## Lens cleaning tissue

- Supplied in packs of 50

| Catalogue No | Dimensions, mm |
|--------------|----------------|
| FB13066      | 80 x 100       |
| FB13067      | 95 x 135       |



## Heat resistant paper

- Heat resistant ceramic paper, approx. 1mm thick
- Dimensions, mm: 500 x 1,000

Catalogue No  
HCR-350-W

 **Bench protection paper, Bench accessories, General Labware, page 385**