


Pipette accessories / Pipettors / Pipettor tips

Pipette stand, polypropylene, vertical



- White, non-corroding stand accommodates 44 pipettes of up to 15mm O.D.
- Bottom tray has small holes for drainage and a separate bush is supplied for height adjustment
- Base diameter, mm: 230
- Rod (dia. x l), mm: 12 x 400

Catalogue No
PMR-541-Y

Pipette jar, polypropylene



- Suitable for storing pipettes, with tight fitting cover
- Autoclavable to 120°C (except polyethylene cover)
- Welded flanged base gives excellent stability
- Can be easily cut down to suit shorter length pipettes
- Can be used as a mini soak jar or canted on side for drainage whilst pipettes in-situ
- Dimensions (dia. x h), mm: 83 x 420 (410)

Catalogue No
PMR-555-010U

Pipette and burettes rinsing set, automatic

- Enables the washing of pipettes up to 600mm long and Mohr burettes to be carried out quickly and safely
- Ideal when radioactive, caustic or toxic materials are in use
- Baskets minimise breakage during washing
- Low water flow rate does not pose a problem
- Extension handle (PMR-569-504K) is for use with basket (PMR-569-502X) and burette rinser (PMR-569-010N)
- Supplied without tubing and clip



Catalogue No	Description	Dia., mm	Height, mm
PMR-565-010N	Burette rinser*	165	990
PMR-567-010B	Pipette rinser**	165	734

Accessories

Catalogue No	Description	Dia., mm	Height, mm
PMR-569-500S	Jar, 165mm x 650mm	165	650
PMR-569-502X	Basket, 130mm x 300mm	130	300
PMR-569-504K	Extension handle	-	-
CNK-641-010F	Tubing clip, for tubing 10 to 14mm diameter, pack of 5	-	-
TWR-670-212T	PVC tubing, 12.5mm x 2.5mm, length 10m	-	-

* Needs 1 basket, 1 jar and 1 extension handle

** Needs 1 basket and 1 jar

Pipettors, dispensing, ECO



- Ideal for schools, colleges and universities
- Repeatability between 0.072% and 0.34%
- Simple to service
- Working accuracy >0.75%
- Rotatable tip ejector for comfort
- Robust and lightweight with controlled spring pressure to cut down hand fatigue
- Continuously variable in their respective capacity ranges
- Coloured for easy differentiation between model sizes



Catalogue No	Volume range, µL	Colour
PMP-705-010N	100 to 1,000	Blue
PMP-705-050B	500 to 5,000	Pink
PMP-705-070S	1,000 to 10,000	Turquoise

Pipettor tips

Catalogue No	Volume, µL	Colour	Pack qty
PMP-675-012X	1 to 200	Yellow	1,000
PMP-675-022L	1,000	Blue	1,000
PMP-675-030M	5,000	Natural	500
PMP-675-040J	10,000	Blue	400