



Pumps, vacuum and pressure, diaphragm, oil-free

- Applications include: vacuum filtration, pressure filtration, analysers; gas sampling analysis and monitoring; circulation and transfer of gases and vapours; flame photometers; vacuum desiccation, small evaporators; pick and place operations; replacement for water jet pumps; vacuum forming
- Range of oil free diaphragm pumps housed in robust metal cases with handle
- Designed specifically for use in educational establishments
- Virtually maintenance free
- Inlet and outlet connectors for 6mm I.D. tubing
- Motors suitable for continuous operation
- Supplied with detachable mains lead

Catalogue No	Flow rate, L/min	Vacuum, mbar	Pressure, bar	Dimensions, ext. (w x d x h), mm	Mass, kg
PXW-610-010X	5	25	0.3	185 x 165 x 155	2.0
PXW-620-010H	6	100	2.4	185 x 165 x 155	1.75
PXW-630-010A	15	100	-	185 x 165 x 155	2.5
PXW-650-010J	39	213	2.0	230 x 160 x 200	4.3

Spares Kits

Catalogue No	Description
PXW-655-500Y	For PXW-610-010X
PXW-655-510F	For PXW-620-010H
PXW-655-520C	For PXW-630-010A
PXW-655-540T	For PXW-650-010J

Pump, vacuum, diaphragm, oil free, MP-Series



- Applications include: filtration; dessication; removal of highly volatile solvents
- PTFE, chemically resistant, oil free vacuum pump for solvent and acid vapours
- Extended warranty, long service intervals, low costs
- No oil to change - ever
- Gas ballast fitted
- Ready to use, with 8mm hose nozzles, mains lead and UK plug

Catalogue No	Flow rate, L/min	Vacuum, mbar
PYF-300-025H	15/30/63	75

**Why not order online?
it's so simple!**



www.griffineducation.co.uk