



Fisherbrand®

## Joint clamps, PTFE

- Chemically inert and totally unaffected by laboratory chemicals
- Temperature resistant, they maintain strength right up to 280°C
- Manufactured from pure PTFE, with totally encapsulated steel springs, they can be adjusted simply by compressing by hand
- Made to fit standard glass joints.
- Supplied in packs of 5



## Connectors, screwthread



- For use with Quickfit® screwthread glass joints
- Comprises of plastic caps, silicone rubber rings and PTFE washers
- Supplied in packs of 10

Catalogue No	Joint size
FB60029	12/-
FB60030	14/-
FB60031	19/-
FB60032	24/-
FB60033	29/-
FB60034	34/-
FB60035	40/-

Catalogue No	Thread size, mm
Plastic caps	
QAK-105-E	13
Rubber rings	
QAK-165-Y	13
QAK-175-B	18
PTFE washers	
QAK-225-F	13
QAK-235-V	18
QAK-237-J	24



## Connectors, screwthread



- The safe and easy way to connect tubing to glass apparatus
- Screws onto glassware fitted with 13mm threads, eg, most Quickfit® condensers
- Supplied in packs of 10

Catalogue No	Description
QAK-350-G	Polypropylene screwthread connector
QAK-353-L	PTFE screwthread connector
QAK-355-W	Spare Viton 'O'-ring



## Connectors, screwthread holder



- For use with connectors, with 13mm screwthread inlets and outlets
- Screws onto the screwthread connector to hold it whilst inserting onto the tubing
- Allows safe and easy removal of the connector from the holes
- Autoclavable
- Supplied in packs of 5

Catalogue No
QAK-357-K