

Gloves, polyethylene, disposable, powder free, non-sterile

- Provide adequate covering for light work
- Can be worn once and thrown away if contaminated
- Ambidextrous
- Supplied in packs of 100 gloves

NOTE: Not to be used with chemicals

Catalogue No	Size
SAR-660-030K	Small
SAR-660-050E	Medium
SAR-660-070V	Large



Gloves, latex, disposable, aloe vera, powder free, **Fisherbrand**® Class I Medical

- Textured surface provides ideal grip for dry and wet working conditions
- Aloe vera coating on the inner surface protects hands from drying, chapping and cracking
- Ideal for prolonged use
- 50mg/g or less of total water extractable protein
- Supplied in packs of 100 gloves

Catalogue No	Size
FB51953	Small
FB51955	Medium
FB51957	Large
FB51959	Extra large



Gloves, vinyl, disposable, lightly powdered, Dura Touch™, PPE Category III

- Clear vinyl seamless gloves with a rolled cuff, AQL 4.0
- Ambidextrous
- Supplied in packs of 100 gloves

Catalogue No	Size
SAR-705-030X	Small
SAR-705-050Y	Medium
SAR-705-070C	Large
SAR-705-090T	Extra large



Gloves, vinyl, disposable, Semperguard, PPE Category I

- Available in lightly powdered and powder free
- Tear resistant, seamless gloves free of latex proteins and accelerators
- Rolled cuff, 0.28mm thickness at the palm and an AQL of 1.5
- Ambidextrous
- Chemical resistance data available on request
- Supplied in packs of 100 gloves

Catalogue No	Size
Lightly powdered	
SAR-746-030D	Small
SAR-746-050U	Medium
SAR-746-070X	Large
SAR-746-090Y	Extra large
Powder free	
SAR-748-030X	Small
SAR-748-050Y	Medium
SAR-748-070C	Large
SAR-748-090T	Extra large