



Stands, burette and funnel

- Hardwood, with interchangeable components
- Bases are 206mm x 110mm (single) and 290mm x 110mm (double)
- Rods are 460mm x 16mm (dia.)
- Burette clamps have lined jaws and funnel supports are slotted
- All are fitted with a plastic screw and rubber feet on base

| Catalogue No | Description |
|--------------|-----------------------|
| STA-200-011D | Burette stand, single |
| STA-200-021A | Burette stand, double |
| STA-200-061L | Funnel stand, single |
| STA-200-071Y | Funnel stand, double |
| STA-200-131Q | Burette/funnel stand |

Spare Parts

| | |
|--------------|------------------------|
| STA-204-010E | Burette clamp, single |
| STA-204-060M | Funnel support, single |
| STA-204-070J | Funnel support, double |



Stands, bases, slotted

- Cast metal base, finished in acid resistant paint designed to support flat objects up to 34mm thick in a vertical position
- Fitted thumb screw
- Dimensions (w x l), mm: 95 x 75
- Mass, kg: 1.4

Catalogue No
STA-395-W



Stands, bases, retort, rectangular

- White acrylic finished steel with hole threaded M10 at the centre of one short side
- Fitted with rubber feet
- Without rod

| Catalogue No | Length, mm | Width, mm |
|--------------|------------|-----------|
| STA-470-030D | 160 | 100 |
| STA-470-050U | 200 | 125 |
| STA-470-070X | 250 | 160 |
| STA-470-090Y | 315 | 200 |



Splints, Gas Analysis, General Labware, page 443
Springs, Properties of Solids, Physics, page 22



Stands, bases, retort, rectangular **Fisherbrand®**

- Pressed steel finished in acid resistant textured acrylic
- Threaded M10
- FB68445 has one hole located centrally on the length and one on the width
- Without rod

| Catalogue No | Description | Length, mm | Width, mm |
|--------------|-------------|------------|-----------|
| FB68445 | Double hole | 315 | 200 |
| FB68401 | Single hole | 160 | 100 |
| FB68403 | Single hole | 200 | 125 |
| FB68405 | Single hole | 250 | 160 |

Bases, retort stand, A shape



- Manufactured from cast iron and painted with a chemically resistant acrylic paint
- The rod mounting hole is dual purpose and will accept both plain rods (secured with supplied thumbscrew) and M10 screwed rods of 12.7mm diameter
- Dimension (w x d x h), mm: 187 x 213 x 67
- Without rod

Catalogue No
STA-630-100U