



Stands, burette clamps


- Burette is held in a three point grip between two formed stops in stainless steel and an eccentric cam, all of which are rubber covered
- Spring loaded cam operated by a lever which is depressed against the spring to open the jaws
- May be used for horizontal and vertical clamping onto rods 10 to 13mm diameter

Catalogue No	Type
STE-800-010J	Single
STE-800-030D	Double

Stands, burette clamp, semi-micro

- Comprises two 13mm plastic coated spring clips mounted 50mm apart on a plate attached to a rod 150 to 8mm diameter, to fit bossheads
- For 25mL or 50mL burettes

Catalogue No
STE-690-070N



Stands, closed rings, with rod

- Plated mild steel with 8mm diameter stem
- Stem length, mm: 140

Catalogue No	Ring dia., mm
STF-520-010J	50
STF-520-031B	70
STF-520-051S	100
STF-520-071M	127
STF-520-091G	150



Stands, closed rings, with bosshead

- With integral bosshead (With plastic headed clamping screw) to accept rods from 8 to 12.5mm diameter

Catalogue No	Dia., mm
STF-550-011J	50
STF-550-031D	75
STF-550-051U	100
STF-550-071X	130



Stands, rings, support plate

- Steel, 100mm diameter, with support rod (dia. x l), 8mm x 125mm

Catalogue No
STF-410-050M



Stands, funnel holders, polypropylene

- Single or double holders for funnels from as small as 25 (with adapter supplied) up to 170mm diameter
- Without the adapter it can be used with powder funnels providing diameter not exceeding 35mm
- Two olives for use with support rods from 8 to 14mm
- Funnel size, mm: 25 to 170
- Rod diameter, mm: 8 to 14

Catalogue No	Places, number of
STF-600-010P	1
STF-600-020M	2