



Syringes, glass, metal luer fitting

Catalogue No	Capacity, mL	Subdivisions, mL
SZR-280-011K	2	0.1
SZR-280-021H	5	0.2
SZR-280-031E	10	0.2
SZR-280-041B	20	1.00
SZR-280-051V	50	2.00



Syringes, luer fitting

- Fully interchangeable

Catalogue No	Capacity, mL	Subdivisions, mL
SZR-260-030U	1	0.05
SZR-260-050X	2	0.1
SZR-260-070Y	5	0.2
SZR-260-090C	10	0.2
SZR-260-110W	20	1.0



Syringes, plastic, luer fitting, disposable, ISO/DIS/7886

- For gases or liquids
- Individually wrapped in sealed envelope
- Polypropylene barrel and plunger, and rubber piston with double contact ring
- Sterilised by ethylene oxide
- Without needles

Catalogue No	Capacity, mL	Type	Pack qty
SZR-205-030J	1 x 0.01	Tuberculin	100
SZR-205-050D	2 x 0.1	-	100
SZR-205-071S	5 x 0.2	-	100
SZR-205-092K	10 x 1.0	Eccentric	100
SZR-205-110L	20 x 1.0	Eccentric	50
SZR-205-130F	30 x 1.0	Eccentric	50
SZR-205-160T	50 x 1.0	Eccentric	25



Syringes, needles, disposable, ISO/DIS/7886

- For SZR-205 and SZR-260 series
- Stainless steel with colour coded polypropylene luer fitting
- Individually packed
- Sterilised by ethylene oxide
- Supplied in packs of 100

Catalogue No	Gauge	Length, mm	Dia., mm	Colour
SZR-275-502X	25	16	0.5	Orange
SZR-275-506G	25	25	0.5	Orange
SZR-275-510P	23	25	0.6	Blue
SZR-275-514H	23	30	0.6	Blue
SZR-275-518W	21	38	0.8	Green
SZR-275-522Y	19	51	1.0	Cream