


Test tube rack, aluminium

Fisherbrand®

Stackable Z-pattern with 5 x 28mm, 6 x 19mm and 8 x 13mm diameter holes, in three rows.



Catalogue No	Description
FB51397	Test tube rack, aluminium



Test tube racks, polypropylene




- Will withstand temperatures up to 120°C
- Indentations on base area to hold test tubes in position


Catalogue No	Shelves	Hole dia., mm	No. of holes	Hole format	Dims. (w x d x h), mm
TEX-360-T	2	19	10	1 x 10	250 x 60 x 74
STK-630-010N	3	13	12	2 x 6	192 x 61 x 81
STK-630-030H	3	21	12	2 x 6	192 x 61 x 81

Test tube racks, polycarbonate

- Slant test tube racks designed for the preparation of two different agar slants in the same rack, end plate allows 5° to 20° slant
- Suitable for use in water baths
- 4 x 10 array with holes keyed with numbers and letters for easy sample identification
- Can be autoclaved occasionally
- Dimensions, mm: 285 x 115



Catalogue No	Tube dia. (max.), mm
STK-635-020D	20



Test tube racks, ResMer™, Unwire™



- Constructed from chemical resistant polyketone, moulded in one piece designed to hold a full load without twisting or bending
- Suitable for use in water baths and freezers
- Moulded letters and numbers on top tier for easy sample identification
- Large end plates allow clear labelling

Catalogue No	Dims. (w x d x h), mm	Tube dia. (max.), mm	Colour	Hole format
STK-770-010A	102 x 200 x 57	13	White	6 x 12
STK-770-012T	102 x 200 x 57	13	Orange	6 x 12
STK-770-014P	102 x 200 x 57	13	Yellow	6 x 12
STK-770-016L	102 x 200 x 57	13	Blue	6 x 12
STK-770-018H	102 x 200 x 57	13	Green	6 x 12
STK-770-020U	102 x 200 x 57	13	Red	6 x 12
STK-770-030R	127 x 248 x 70	16	White	6 x 12
STK-770-032N	127 x 248 x 70	16	Orange	6 x 12
STK-770-034J	248 x 127 x 70	16	Yellow	6 x 12
STK-770-036F	248 x 127 x 70	16	Blue	6 x 12
STK-770-038B	248 x 127 x 70	16	Green	6 x 12
STK-770-040X	248 x 127 x 70	16	Red	6 x 12
STK-770-050L	250 x 102 x 83	20	White	4 x 10
STK-770-052H	250 x 102 x 83	20	Orange	4 x 10
STK-770-054D	250 x 102 x 83	20	Yellow	4 x 10
STK-770-056W	250 x 102 x 83	20	Blue	4 x 10
STK-770-058S	250 x 102 x 83	20	Green	4 x 10
STK-770-060Y	250 x 102 x 83	20	Red	4 x 10
STK-770-070F	300 x 121 x 92	25	White	4 x 10
STK-770-072B	300 x 121 x 92	25	Orange	4 x 10
STK-770-074U	121 x 300 x 92	25	Yellow	4 x 10
STK-770-076Q	121 x 300 x 92	25	Blue	4 x 10
STK-770-078M	121 x 300 x 92	25	Green	4 x 10
STK-770-080C	121 x 300 x 92	25	Red	4 x 10
STK-770-090W	108 x 283 x 83	30	White	3 x 8
STK-770-092S	108 x 283 x 83	30	Orange	3 x 8
STK-770-094X	108 x 283 x 83	30	Yellow	3 x 8
STK-770-096K	108 x 283 x 83	30	Blue	3 x 8
STK-770-098G	108 x 283 x 83	30	Green	3 x 8
STK-770-110T	108 x 283 x 83	30	Red	3 x 8