



Timer, count down, jumbo display, digital

- Counts down/up from 99min 59sec with alarm
- 27mm and 20mm digit height
- Back clip with stand, magnet and hook
- Dimensions, mm: 90 x 65
- Supplied with battery

Catalogue No
TKP-150-K

Spare Part	Description
Catalogue No BMT-680-040W	Alkaline batteries 1.5V (AAA), pack of 4



Timer, count down triple timer, digital

- Counts down from 20hr with alarm
- Memory timer
- 1/100s stopwatch
- Clock
- 10mm digit height
- Back clip with stand and magnet
- Dimensions, mm: 60mm x 70mm

Catalogue No
TKP-394-L

Spare Part	Description
Catalogue No BMD-110-021W	Alkaline battery 1.5V (LR44)

Timer, count down/up, digital

- Counts down/up from/to 10hr housed in a robust plastic case with magnetic clip
- Large easy to read LCD
- 12/24hr clock
- Supplied with battery

Catalogue No
TKP-340-B



Timer, count down, digital

- Counts down from 99min 59s with alarm housed in a robust case with magnet
- Large easy to read LCD
- Dimensions (w x d x h), mm: 46 x 12 x 88
- Mass, g: 30



Catalogue No
TKP-440-Q

Timer, controller, 24 hour, digital

- Easily programmable slimline digital timer controller providing accurate to the minute ON/OFF control for electrical appliances
- Easy to view display allows up to four ON/OFF settings per day and incorporate a self cancelling override for instant switching
- For 230V, 50Hz single phase supplies



Catalogue No
TKR-200-R