



## Gauges

- Vernier calliper, general purpose with thumb wheel fine adjustment
- Can be used for internal, external and depth measurement
- Metric, imperial dual graduations


Catalogue No	Description
TOH-661-010N	Length 115mm
TOH-661-020K	Length 150mm
TOH-661-030H	Length 150mm, with plastic body
TOH-700-K	Micrometer screw gauge, 0 to 25mm x 0.01mm



## Glass cutters

Catalogue No	Description
TOL-130-010R	① Cutter, glass tubing, suitable for soda and borosilicate glass tubing up to 40mm dia.
TOL-140-N	Cutter, glass tubing, manufactured in nickel plated zinc alloy and accepts tubing up to 30mm dia.
TOL-135-010K	Diamond, glass cutting, in metal mount ② with wooden handle
GTB-490-R	Glass knife, super-hard, wooden handle, ③ blade, 65mm x 13mm x 1.6mm

Spare Parts	
TOL-144-010J	Wheels for TOL-140-N, pack of 5



## Gauges, digital

**Fisherbrand®**

- Measures I.D, O.D, and depth
- Measures in inches or mm
- Moisture and dust resistant steel
- Jumbo digit display
- Supplied with an ISO 17025, A2LA Traceable® Certificate to US NIST (National Institute of Standards and Technology)

Catalogue No	Description
FB70250	150mm to 240mm
FB70252	200mm to 320mm



## Hammers

Catalogue No	Description
TOH-743-K	Claw, wooden handle, 16oz
TOH-761-Y	Steel claw, 8oz
TOH-763-T	Steel claw, 16oz



## Hole punch

Catalogue No	Description
TOH-640-N	Hole punch, eyeletter
TOH-644-010J	Replacement eyelets