



Tubing, soda-lime **Fisherbrand** glass

- Tolerance $\pm 0.2\text{mm}$: wall thickness 0.7mm
- In lengths of approximately 0.5m
- Supplied in packs of 50

Catalogue No	Dia. (ext.), mm
FB51460	4
FB51461	5
FB51462	6
FB51463	7
FB51464	8
FB51465	10

Tubing, borosilicate glass **Fisherbrand**

Catalogue No	Ext. dia., mm	Bore size, mm
Medium wall, pack of 20		
FB51467	6	3
FB51468	7	4
FB51469	8	5
FB51470	10	7
Capillary bore, pack of 10		
FB51471	5	0.6
FB51472	6	0.8
FB51473	7	1.2




Tubing connectors, polypropylene, T-shape



- Conforms to food and agricultural regulations
- Dimensions are nominal

Catalogue No	O.D., mm	Bore size, mm
ADF-715-030L	5.0	2.8
ADF-715-050F	7.0	4.4
ADF-715-070W	9.0	6.3
ADF-715-090Q	11.0	8.2
ADF-715-110N	13.0	10.0
ADF-715-010R	3.0	1.6
ADF-715-130H	15.0	12.1



Tubing clips

Catalogue No For tubing up to Pack
to dia., mm qty

Ratchet, polypropylene ①
Suitable for use with flexible tubing up to 10mm diameter

CNK-341-J	10	10
-----------	----	----

Open sided ②
Plated metal screw clip with broad pressure plate for blood transfusion use

CNK-381-E	13	10
-----------	----	----

Hoffman ③ ④
Plated nickel, with clamping screw hinges

CNK-410-011A	15	10
CNK-410-031R	20	10
CNK-410-051L	25	10
CNK-410-071F	32	10
CNK-410-092U	40	10

Hoffman, heavy duty ⑤
Nickel plated and extra strong, for heavy wall flexible tubing, with hinged bottom bar

CNK-440-H	35	each
-----------	----	------

Screw, tubing ⑥
Moulded in 'Delrin' acetal homopolymer, the hinged side allows it to be clamped around tubing already connected to apparatus

CNK-481-T	14	10
-----------	----	----

Dialysis ⑦
Plastic, for closing the end of dialysis visking tubing up to 38mm inflated diameter

CNK-490-S	-	10
-----------	---	----

Adjustable flow ⑧
Made from Delrin® polyacetal plastic, with easy to operate roller control for liquids and gases, non-slip serrated wheel and sloped floor allows infinitely variable flow regulation, heat resistant up to 150°C

CNK-500-032X	4.5	5
CNK-500-072C	10	5
CNK-500-092T	14	5



Tubing, clips, spring

- Made from chemical and stain resistant acetal copolymer and stainless steel for compressing flexible tubing

Catalogue No
CNK-700-010Y



Tubing holder



- Holds flexible tubing steady
- Clips onto beakers up to 8mm thick
- Suitable for tubing up to 9.5mm diameter
- Supplied in packs of 6

Catalogue No
CNK-320-001E



Tubing connectors, polypropylene, tapered



- For straight connection of different diameter tubing

Catalogue No	Dia. (ext.), mm	Dia. (int.), mm	Pack qty
ADF-642-030Q	3 to 5	1.9	20
ADF-642-050K	4 to 8	2.8	20
ADF-642-070E	6 to 10	3.9	20
ADF-642-090V	8 to 12	6.0	20
ADF-642-110S	10 - 14	7.9	10
ADF-642-130M	12 - 16	10.0	10