




Test tubes, borosilicate glass, medium wall, rimmed

page
544

Catalogue No	Length, mm	Dia. (ext.), mm	Wall thickness, mm	Pack qty
TES-610-030C	75	10	1.0	100
TES-610-050T	75	12	1.0	100
TES-610-070N	100	12	1.0	100
TES-610-090H	100	16	1.2	100
TES-610-110E	125	16	1.2	100
TES-610-130V	150	16	1.2	100
TES-610-170J	150	18	1.2	100
TES-610-210A	150	24	1.2	50




Test tube rack, hardwood

- Single row, hardwood top and bottom members with moulded polypropylene supports which clip into position
- Occupies only about one-third of the space required by a conventional rack when dismantled for storage
- Supplied in packs of 5
- No. of holes: 4 x 22mm/2 x 32

Catalogue No
TEX-304-030F ①

page
547

Test tube rack, aluminium **Fisherbrand®**




Stackable Z-pattern with 5 x 28mm, 6 x 19mm and 8 x 13mm diameter holes, in three rows

Catalogue No
FB51397

test tubes

page
546

Test tube racks



Catalogue No
STK-750-010X ①

Dimensions (w x d x h), mm
200 x 78 x 60

Hole format
2 x 6

Nylon coated steel, tubes ②

Catalogue No
STK-640-030A
STK-640-050R
STK-640-130T

Places,
no. of

Tube dia.
(max.), mm
10 to 14
13 to 18
16 to 25

Tube length
(min.), mm
50
100
100

page
546

Test tube holders



Catalogue No Description

Wood with steel spring ①
TEW-650-M Wood with steel spring, for tubes up to 19mm (dia.), pack of 10

Plated spring steel ②
TEW-710-J Plated spring steel, for tubes up to 30mm (dia.), pack of 10

Call: 01509 233344

Fax: 01509 231893