


### Inoculating spreaders, polystyrene, sterile

**S** 03



High impact polystyrene disposable bacterial spreaders.

Wedge shaped to reduce tendency to flip or traditional 'L' shape. Packs of 500 contain re-sealable inner bags of either 20 or 5 sterile spreaders.

Wedge shaped

Catalogue No	Inner pack qty	Pack qty
<b>LPS-141-011R</b>	20	500
<b>LPS-141-021X</b>	5	500

L-shaped

Catalogue No	Inner pack qty	Pack qty
<b>LPS-140-041X</b>	5	500
<b>LPS-140-031R</b>	20	500

### Inoculating spreaders, bacterial, cell

**scienceware**

**a** 04



Stainless steel cell spreaders with a smooth triangular end and a straight spreading bar.

Autoclavable and chemically sterilisable. Both sizes are 203mm long. Handle diameter 4.8mm.

Catalogue No	Description
<b>MGF-200-010J</b>	Small cell spreader for 60mm plates has a spreading bar 16mm wide
<b>MGF-200-030D</b>	Large cell spreader for 90mm plates has a spreading bar 25mm wide

### Inoculating spreaders, polystyrene, sterile

**Fisherbrand**

**S** **NEW** 46



High impact disposable spreader for spreading over the surface of solid media for microbiological applications.

- Packed in tamper evident zipper bags
- L shape spreaders have curved ends
- T shape spreader for more precise inoculating

Catalogue No	Description	Inner pack qty	Pack qty
<b>FB58782</b>	L shape, 148mm	1	400
<b>FB58781</b>	L shape, 148mm	5	400
<b>FB58780</b>	L shape, 148mm	25	1,000
<b>FB58836</b>	T shape, 150mm	25	3,500