

Pipettes, serological, glass, non-sterile



02

- Non-sterile, packed in plastic sleeves
- Graduated glass pipette with black graduations

Accuracy is within $\pm 2\%$ of volume, which is adjusted for blowout. 10mL size has polypropylene top stem of reduced diameter, this enables finer control to be achieved when hand pipetting and allows superior plugging. **PMK-100-020A** and **PMK-100-030U** supplied with PVT teat and adapter.



Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Length, mm	Jet	Inner pack qty	Pack qty
PMK-100-020A	D804	1	0.01	230	Short	10	500
PMK-100-030U	D805	1	0.10	230	Short	10	500
PMK-100-040R	D802	1	0.01	270	Drawn	10	500
PMK-100-050X	D806	2	0.02	260	Short	10	500
PMK-100-060L	D803	2	0.02	270	Drawn	10	500
PMK-100-070Y	D807	5	0.05	310	Drawn	10	500
PMK-100-080F	D808	10	0.10	310	Drawn	10	250

Pipettes, serological, sterile



02

- Graduated glass pipettes

Catalogue No	Alt. No	Capacity, mL	Graduations, mL	Length, mm	Jet	Inner pack qty	Pack qty
PMK-200-040J	D804/PS	1	0.01	230	Short	10	500
PMK-200-041H	D805/PS	1	0.10	230	Short	10	500
PMK-200-042F	D802/PS	1	0.01	270	Drawn	10	500
PMK-200-050G	D806/PS	2	0.02	260	Short	10	500
PMK-200-052C	D803/PS	2	0.02	270	Drawn	10	500
PMK-200-070A	D807/PS	5	0.05	310	Drawn	5	500
PMK-200-080U	D808/PS	10	0.10	310	Drawn	5	250



Pipette filler

Fisherbrand®


46

- Use with all glass and plastic pipettes
- Low battery indicator
- Adapters for 5mL and 10mL pipettor tips
- Support allows unit to be rested with pipette attached

A lightweight, cordless, rechargeable pipetting aid which can be used with all plastic and glass pipettes between 1mL and 100mL.

Linear speed control allows the aspiration and dispensing speed to be adjusted easily and precisely to suit either large or small volume pipettes. The speed can also be fine tuned by applying varying finger pressure to the operating buttons. The fold out stand can be used to hold the unit and pipette when not in use.

Technical Specification

Filter	Hydrophobic 0.45µm
Battery	3 x AA 1.2V/270mAh, NiCd.
Battery life, hr	7 (Approx.)
Electrical supply	a.c. charger, 9V d.c. (input 230V, 50/60Hz). Charging time: 14hr when fully discharged

Catalogue No	
FB69000	

Spare parts

Catalogue No	Description
PMP-470-010S	Proline tip 5000 plus, in box, 150mm, 100µL to 5,000µL, pack of 100
FB69014	Replacement filter, non-sterile, pack of 5
FB69005	Replacement filter, sterile
FB69006	5mL nose cone
FB69007	10mL nose cone
FB69008	25mL nose cone
FB69009	5mL silicone adapter
FB69010	10mL silicone adapter
FB69011	25mL silicone adapter
FB69012	Battery pack

