

Cell/Tissue Culture and Cell Transfection

Apparatus and consumables - TC-treated tubes

Tubes, cell culture, polystyrene, Thermo Scientific Nunc®, Nuncion™ Δ



76

- Certified surface treatment for optimal cell attachment and growth
- Flat sided tube for adherent cells
- Flat side allows easy microscopy, a cover slip can be used in the tube
- Medium retained in the flat sided tube in horizontal position
- Excellent optical quality
- Certified surface treatment for optimal cell attachment and growth
- Sterile



Catalogue No	Alt. No	Shape/ configuration	Cap style	Suggested working volume, mL	Culture area, cm ²	Dimensions, external, mm	Inner pack qty	Pack qty
TKT-150-030H	146183	Round	Screw cap, PE	7	—	100 x 14	100	600
TKT-150-010N	145470	Round	Push-on cap, PE	7	—	100 x 15	100	600
TKT-150-050B	156758	Flat	Screw cap, PE	3	5.5	110 x 16	75	450

Tubes, cell culture, polystyrene, TC-treated

CORNING



66

- Threaded plug seal cap to prevent leakage
- Optically clear
- 15mL capacity
- Dimensions: 16 x 125mm
- Sterilised by gamma radiation
- Certified non-pyrogenic

Catalogue No	Alt. No	Treated, [Y/ N]	Cap style	Pack type	Inner pack qty	Pack qty
TKV-180-056R	430172	Yes	Screw top	Rack	50	500
TKV-180-110X	430157	No	Screw top	Bulk	25	500



Tubes, cell culture, plastic, TC-treated



02

- Suitable for suspension culture procedures
- Sterile bulk packs
- Push on caps
- Dual position cap, loose for gas exchange and push-down for tight seal
- With large white marking spot for easy labelling
- Printed graduations

Polystyrene

Catalogue No	Alt. No	Dimensions [dia. x h], mm	Packaging type	Inner pack qty	Pack qty
TKV-182-011E	75CTS1	12 x 75	Individually wrapped	1	500
TKV-182-015T	75CTS	12 x 75	Bulk packed	25	500
TKV-182-051P	100CTS1	17 x 100	Individually wrapped	1	500
TKV-182-055H	100CTS	17 x 100	Bulk packed	25	500

Polypropylene

Catalogue No	Alt. No	Dimensions [dia. x h], mm	Packaging type	Inner pack qty	Pack qty
TKV-182-031V	75CTP1	12 x 75	Individually wrapped	1	500
TKV-182-035N	75CTP	12 x 75	Bulk packed	25	500
TKV-182-071J	100CTP1	17 x 100	Individually wrapped	1	500
TKV-182-081G	100CTP	17 x 100	Bulk packed	25	500

