



Please contact Fisher Scientific to discuss your EU compatible product requirements.

Accessories for CryoBank™ and Bank-IT™



76

Handheld electric screwdriver

- The cap is retained by the socket, allowing single handed manipulation

Manual screwdriver

- For manual unscrewing and fastening of caps

Bench top rack

- Holds plate whilst vials are manipulated, filed and emptied

Tube selection

- For quick picking of multiple vials, avoids thawing of remaining vials

Technical Specification - Specific

	CRY-369-644S PJB-510-020P	CRY-369-658H PJB-510-010S
Electrical supply	220/240V 50/60Hz	240V 50Hz

Catalogue No	Alt. No	Description	Pack qty
CBT-300-200Y	374009	Press out tool, Bank-IT™	1
CBT-300-210F	374011	Retention plate, Bank-IT™ (rack holder)	1
CRY-369-640D	369640	Cryobox 13 x 13 divider cardboard for CryoBank™ and Bank-IT™ vials	48
CRY-369-643U	369643	Screwdriver small hand held for CryoBank™ and Bank-IT™ vials black	1
CRY-369-644S	369644	Screwdriver electric european plug 220-240V 50/60Hz for CryoBank™ and Bank-IT™ vials	1
CRY-369-658H	369658	Screwdriver electric UK plug 240V 50Hz for CryoBank™ and Bank-IT™ vials	1
PJB-510-010S	330042	2D reader Nunc rack scanner, EU plug	1
PJB-510-020P	330043	2D reader Nunc rack scanner, UK plug	1



Tube racks, cluster, grid-referenced and twist-lock, polypropylene, Thermo Scientific Matrix SmartScan™



04

For sample storage and sample tracking. The SmartScan™ system is designed for high throughput, individual sample identification.

- Data Matrix™ and human readable code on each tube with over 3.6 x 10⁹ unique code combinations
- 96 individual tubes
- Encapsulated code for protection against damage
- Temperature range: -180°C to 120°C
- Heavy duty rack with lockable microplate format lid SBS footprint
- Each tube is fully supported to withstand pressure applied by robotic systems
- Codes can be scanned whilst tubes are racked
- Alphanumeric grid-referencing moulded into the rack lid
- Tubes are available racked or loose in bags
- Available unsealed or sealed with MULTISIP™ septum plugs
- Barcoding available
- Contamination free
- Cleanroom produced