

### Microplate, storage, 384 well, 120µL, Thermo Scientific ABgene

**Thermo**  
SCIENTIFIC

03

- Maximum well volume of 120µL
- Pyramidal well bottom for ease of sample retrieval
- Discrete raised rims around each well for effective sealing and to reduce the risk of cross contamination
- Additional raised corner rims to protect edges of seals

Catalogue No	Alt. No	Colour	Pack qty
<b>TUL-050-100N</b>	AB-0781	Natural	50
<b>TUL-050-101L</b>	AB-0781/K	Black	50
<b>TUL-050-108U</b>	AB-0781/Y	Yellow	50



### Microplates, storage, 384 well, polypropylene

CORNING

**S**

66

384 well plates and blocks

- Sterile
- Round bottoms for optimal handling
- Rounded square well geometry reduces wicking
- Resistant to many common organic solvents (e.g., DMSO, ethanol, methanol)
- Certified DNase and RNase free
- Both plates and blocks are automation friendly

384 well plates

Catalogue No	Alt. No	Colour	Well shape	Sterile, Yes/No	Well volume, µL	Inner pack qty	Pack qty
<b>DPS-164-010X</b>	3656	Clear	Round	Yes	95	25	100
<b>DPS-164-020L</b>	3657	Clear	Round	No	95	25	100
<b>DPS-164-050C</b>	3658	Black	Round	No	95	25	100

384 well blocks

Catalogue No	Alt. No	Colour	Well shape	Sterile, Yes/No	Well volume, µL	Inner pack qty	Pack qty
<b>DPS-164-030Y</b>	3964	Clear	Round	Yes	180	25	25
<b>DPS-164-040F</b>	3965	Clear	Round	No	180	100	100
<b>DPS-164-060W</b>	3342	Clear	V bottom	Yes	240	50	50
<b>DPS-164-070T</b>	3347	Clear	V bottom	No	240	50	50



### Microplates, storage, 384 deep well, polypropylene, Thermo Scientific ABgene

**Thermo**  
SCIENTIFIC

03

Sample storage

- 384 square wells, maximum well volume 300µL
- Pyramidal well bottom for ease of sample retrieval
- Discrete raised rims around each well for effective sealing and to reduce risk of cross contamination
- Additional raised corner rims to protect edges of seals
- Sealing options: heat sealing and adhesive sealing
- Alphanumeric grid-referencing
- Cleanroom produced
- SBS footprint
- Barcoding available
- Contamination free

Catalogue No	Alt. No	Description	Pack qty
<b>TUL-050-120H</b>	AB-0862	384 deep well storage plate	50
<b>TUL-050-123B</b>	AB-0862/a	384 deep well storage plate in amber	50

