



Microplates, PCR, polycarbonate, non-sterile, 96 well, Thermowell®

CORNING

66

- Thermally stable up to 135°C, preventing warping and ensuring a precise fit in the heat block throughout the entire cycle program
- Plate is designed to precisely fit into a specific thermal cycler
- Provides optimum heat transfer
- Certified RNase/DNase free
- Plates are available to fit all current 96 well heat blocks on the market
- Thermowell® has a ridge inside the wells for sample loading and for sequencing reactions
- For oil-free operation, plates must be used with sealing tape or a sealing mat and used in conjunction with thermal cyclers that have heated lids
- Individually wrapped

Catalogue No	Alt. No	Model No	Well volume, µL	Pack qty
MPA-120-010A	6509	P	200	25
MPA-120-020U	6510	C	200	25
MPA-120-030R	6511	M	200	25
MPA-120-040X	6512	T	170	25
MPA-120-060Y	6514	H	300	25



Microplate, PCR, polycarbonate, non-sterile, 96 well, Thermowell®, DNA-BIND™

CORNING

66

Can be used in several solid phase applications, such as DNA synthesis, conventional PCR* reactions, and standard DNA hybridisation assays.

- Individually wrapped
- Fits most commonly used thermal cyclers
- For oil-free operation, plates must be used with sealing tape or a sealing mat and used in conjunction with thermo cyclers with heated lids
- With DNA-BIND™ surface
- With clear conical bottom

Catalogue No	Alt. No	Model No	Well volume, µL	Pack qty
MPA-122-010L	6573	M - microplate, clear	200	25

Microplate accessories, PCR*, lid/plate holder, Thermowell™/Stripwell™

CORNING

66

Lid: Thermally stable up to 135°C and can be used on Thermowell™ 96 well polycarbonate plates as an additional protection to oil overlays.

Plate holder: 96 well egg-crate style holder for 1 x 8 Stripwell™ strips. Holds twelve. Allows each individual well to be positioned back into place once broken apart.

*Polymerase Chain Reaction (PCR) is a process covered by patents owned by Hoffman-La Roche

Catalogue No	Alt. No	Description	Pack qty	Case qty
MPA-124-010W	6520	Fits 6510, 6511 and 6514 microplates	1	25
MPA-124-070E	2572	Plate holder	5	20

For a full range of pipettes and accessories, please see the Fisher Scientific Laboratory Catalogue or visit our website.

