



0.2mL gravity flow columns



1mL gravity flow columns



5mL and 10mL gravity flow columns



1.7mL HPLC columns

Continuation

Strep-tag® protein purification system



LE

Buffers and reagents for Strep-tag®

Catalogue No	Alt. No	Description	Quantity
IB20204050	2-0204-050	Avidin (to block biotin)	50mg
IB20401002	2-0401-002	Anhydrotetracycline (inducer for tet promoter)	25mg
IB20501002	2-0501-002	Biotin blocking buffer	2mL
IB21000002	2-1000-002	D-Desthiobiotin	1g
IB21000025	2-1000-025	D-Desthiobiotin (10x buffer E)	25mL
IB21002001	2-1002-001	Strep-tag® protein purification buffer set	100mL 10x Buffer W 25mL 10x Buffer E 100mL 10x Buffer R
IB21002100	2-1002-100	Strep-tag® regeneration buffer with HABA (10x buffer R)	100mL
IB21003100	2-1003-100	Strep-tag® washing buffer (10x buffer W)	100mL
IB21017050	2-1017-050	IBA-lyse bacterial lysis buffer	50mL
IB21018002	2-1018-002	Strep-tag® II peptide	1.8mg
IB21019025	2-1019-025	Buffer BE (10x, with biotin)	25mL
IB21020100	2-1020-100	Buffer L (5x, for mammalia)	100mL

Column, cartridge and plate formats for Strep-Tactin® resins

(For applied volume, adjust protein extract volume according to column binding capacity)

Catalogue No	Alt. No	Type	Description	Applications	Volume
IB21202550	2-1202-550	Strep-Tactin® Sepharose®	0.2mL gravity flow columns, x 5	Gravity flow	0.1 to 2mL (applied)
IB21207550	2-1207-550	Strep-Tactin® Superflow®	0.2mL gravity flow columns, x 5	Gravity flow	0.1 to 2mL (applied)
IB21209550	2-1209-550	Strep-Tactin® Superflow® high capacity	0.2mL gravity flow columns, x 5	Gravity flow	0.1 to 2mL (applied)
IB21506550	2-1506-550	Strep-Tactin® MacroPrep®	0.2mL gravity flow columns, x 5	Gravity flow	0.1 to 2mL (applied)
IB21202001	2-1202-001	Strep-Tactin® Sepharose®	1mL gravity flow column	Gravity flow	0.5 to 10mL (applied)
IB21207001	2-1207-001	Strep-Tactin® Superflow®	1mL gravity flow column	Gravity flow	0.5 to 10mL (applied)
IB21209001	2-1209-001	Strep-Tactin® Superflow® high capacity	1mL gravity flow column	Gravity flow	0.5 to 10mL (applied)
IB21506001	2-1506-001	Strep-Tactin® MacroPrep®	1mL gravity flow column	Gravity flow	0.5 to 10mL (applied)
IB21202051	2-1202-051	Strep-Tactin® Sepharose®	5mL gravity flow column	Gravity flow	2.5 to 50mL (applied)
IB21207051	2-1207-051	Strep-Tactin® Superflow®	5mL gravity flow column	Gravity flow	2.5 to 50mL (applied)
IB21209101	2-1209-101	Strep-Tactin® Superflow® high capacity	5mL gravity flow column	Gravity flow	2.5 to 50mL (applied)
IB21506051	2-1506-051	Strep-Tactin® MacroPrep®	5mL gravity flow column	Gravity flow	2.5 to 50mL (applied)
IB21202101	2-1202-101	Strep-Tactin® Sepharose®	10mL gravity flow column	Gravity flow	5 to 100mL (applied)
IB21207101	2-1207-101	Strep-Tactin® Superflow®	10mL gravity flow column	Gravity flow	5 to 100mL (applied)
IB21506101	2-1506-101	Strep-Tactin® MacroPrep®	10mL gravity flow column	Gravity flow	5 to 100mL (applied)
IB21203017	2-1203-017	Strep-Tactin® POROS 20	1.7mL HPLC column	HPLC	0.85 to 17mL (applied)
IB21205017	2-1205-017	Strep-Tactin® POROS 50	1.7mL HPLC column	HPLC	0.85 to 17mL (applied)
IB21231001	2-1231-001	Strep-Tactin® Superflow®	1mL H-PR cartridge	FPLC/HPLC	0.5 to 10mL (applied)