



### Microscope slide incubation box, Slideshow

04

Can be used in refrigerators and incubators, -10°C to 45°C

- Neat stackable design to maximise available space
- Range of sizes - 10, 20, 30 slides
- Clear or black lids for fluorescence work
- Wash down and drain off facility for reagents

It is an essential aid in preventing loss or spoilage of histology sections.

The Slideshow is a purpose built, vacuum formed humid incubation box for those who handle slides on a regular basis.

Technical Specification - Specific

	Slideshow10	Slideshow20	Slideshow30
Dimensions [w x d x h], mm	385 x 160 x 50	385 x 245 x 50	385 x 330 x 50
Mass, kg	0.8	1.5	2.1

Technical Specification - General

Spacing [max], mm	35 (slide rail)
Temperature, °C	45 (max.) -10 (min.)

Catalogue No	Description
<b>INH-110-010H</b>	Slideshow10, with clear lid - for up to 10 slides
<b>INH-110-015U</b>	Slideshow10, with black lid - for up to 10 slides
<b>INH-110-020E</b>	Slideshow20, with clear lid - for up to 20 slides
<b>INH-110-025R</b>	Slideshow20, with black lid - for up to 20 slides
<b>INH-110-030B</b>	Slideshow30, with clear lid - for up to 30 slides
<b>INH-110-035X</b>	Slideshow30, with black lid - for up to 30 slides

Spare parts

Catalogue No	Description
<b>INH-110-500N</b>	Clear lid for Slideshow10
<b>INH-110-510K</b>	Black lid for Slideshow10
<b>INH-110-520H</b>	Clear lid for Slideshow20
<b>INH-110-525U</b>	Black lid for Slideshow20
<b>INH-110-530E</b>	Clear lid for Slideshow30
<b>INH-110-535R</b>	Black lid for Slideshow30



### Microscope slide storage tray, polystyrene and ABS

02

**Kartell**

Ideal for storage of large quantities of 76mm x 26mm microscope slides. Each cavity is numbered. Batch data can be written on either the inside or outside of the end plates.

It is strongly recommended that the covers **MNK-384-500H** be used to prevent the ingress of dust and other foreign bodies. Holds 100 slides.

Catalogue No	Alt. No	Width, mm	Depth, mm	Height, mm	Pack qty
<b>MNK-380-070H</b>	00921	38	360	100	2

Accessory dust cover

Catalogue No	Alt. No	Width, mm	Depth, mm	Height, mm	Pack qty
<b>MNK-384-500H</b>	00920	35	345	40	2



### Microscope slide cabinet

46

**Fisherbrand®**

Manufactured from steel, the slide cabinet holds 100 slide holders, a total of 400 slides.

The slide holders (**FB59511**) are manufactured from sturdy cardboard with a protective plastic frame. Fully compatible with all makes and thicknesses of slides, each holder accepts four slides of 76mm x 26mm.

Catalogue No	Dimensions, external [w x d x h], mm
<b>FB59512</b>	128 x 165 x 413

Accessories

Catalogue No	Description	Pack qty
<b>FB59511</b>	Microscope slide holder	12